










DACLON

NYLON
BLACK

		Metric EP	USP	Länge	Bestellnummer
 <p>Rundkörpernadel 2 Nadeln</p> <p>1/4 Kreis</p>	<p>2 x VR 8 mm 220 µm</p> 	0,4	8/0	45 cm	ON81
	<p>2 x HSP 6,2 mm 150 µm</p> 	0,2	10/0	30 cm	ON102
 <p>Spatulanadel 2 Nadeln</p> <p>1/2 Kreis</p>	<p>2 x HSP 6,4 mm 220 µm</p> 	0,2	10/0	30 cm	ON106
	<p>2 x DSP 6,2 mm 150 µm</p> 	0,2	10/0	30 cm	ON101
 <p>Spatulanadel 2 Nadeln</p> <p>3/8 Kreis</p>		0,3	9/0	30 cm	ON95
	<p>2 x DSP 6,4 mm 220 µm</p> 	0,2	10/0	30 cm	ON105
	<p>2 x VSP 8 mm 220 µm</p> 	1	5/0	45 cm	ON51

Augenchirurgie