

DACLON NYLON

Type	Monofilament
Composition	Extrusion of polyamide 6.0 or 6.6
Coating	None
Colour	Blue Black
Tissue reaction	Minimal
Absorption	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
Presentation	Box with 12 sutures

DACLON

NYLON

CHARACTERISTICS

High tensile strength, non absorbable

Synthetic

Monofilament with microscopic uniform diameter

Pliability and softness

Hermetically sealed packing

BENEFITS

Excellent and extended tissue support

Minimal tissue reaction

Soft passage through the tissues without sawing, tissue drag and trauma

Absence of capillarity (no support of bacterial growth)

Excellent knot tie down

Good handling properties

Easy and secure knotting

Guaranteed seal and product sterility

Indications General soft tissue closing and/or ligation

Daclon Blue especially in Skin closure
General surgery
Plastic surgery
Gastrointestinal surgery
Gynaecology
Obstetrics
Orthopaedics

Daclon Black especially in Plastic surgery
Ophthalmic surgery