

		Metric EP	USP	Lenght	Order number
 <p>Side cutting spatula double armed</p> <p>1/2 Circle</p>	<p>2 x HSP 8 mm 200 µm</p> 	0,5	7/0	45 cm	OSU07
 <p>Side cutting spatula double armed</p> <p>3/8 Circle</p>	<p>2 x DSP 6 mm 200 µm</p> 	0,4	8/0	45 cm	OSU08
 <p>Side cutting spatula double armed</p> <p>1/4 Circle</p>	<p>2 x VSP 8 mm 300 µm 350 µm</p> 	0,7 1	6/0 5/0	45 cm 45 cm	OSU06 OSU05