

SUPRAMID

BLACK

SUPRAMID
BLACK

Type	Twisted and coated multifilament
Composition	Polyamide 6.6 strands enclosed in a polyamide 6 sheet
Colour	Black
Tissue reaction	Minimal
Absorption	Non absorbable, gradually encapsulated by connective tissue
Presentation	Box with 12 sutures

The logo for SMI, featuring the letters 'smi' in a stylized, rounded, lowercase font. The 's' is white with a blue outline, and the 'm' and 'i' are solid blue. The 'i' has a small circle above it.

www.sutures.be

SUPRAMID

BLACK

CHARACTERISTICS

High tensile strength, non absorbable

Synthetic

Twisted structure

Multifilament in a sheath

Pliability and softness

Hermetically sealed packing

Indications General soft tissue closing
and/or ligation
especially in Skin closure

BENEFITS

Excellent and permanent tissue support

Minimal tissue reaction

Excellent handling properties

High flexibility

High tensile strength

Excellent knot security

Soft passage through the tissues without sawing, tissue drag
and trauma

Good knot tie down / adjustability

Reduced capillary action

Good handling properties

Easy and secure knotting

Guaranteed seal and product sterility