



Sutures

in cassette packs



*Sutures in
cassette packs*

SURGICRYL®

POLYGLYCOLIC ACID

Braided multifilament
Absorbable



Metric EP	USP	Lenght	Order number
1,5	4/0	15 m	11151
2	3/0	15 m	11201
2	3/0	25 m	11202
3	2/0	15 m	11301
3	2/0	25 m	11302
3	2/0	50 m	11305
3,5	0	15 m	11351
3,5	0	25 m	11352
3,5	0	50 m	11355
4	1	15 m	11401
4	1	25 m	11402
4	1	50 m	11405
5	2	15 m	11501
5	2	25 m	11502
5	2	50 m	11505
6	3&4	15 m	11601
6	3&4	25 m	11602
6	3&4	50 m	11605
7	5	15 m	11701
7	5	25 m	11702
7	5	50 m	11705
8	6	15 m	11801
8	6	25 m	11802
8	6	50 m	11805

SURGICRYL®

MONOFILAMENT

POLYDIOXANONE

Monofilament
Absorbable



Metric EP	USP	Lenght	Order number
1,5	4/0	15 m	13151
2	3/0	15 m	13201
2	3/0	25 m	13202
3	2/0	15 m	13301
3	2/0	25 m	13302
3,5	0	15 m	13351
3,5	0	25 m	13352
4	1	15 m	13401
4	1	25 m	13402
4	1	50 m	13405
5	2	15 m	13501
5	2	25 m	13502
5	2	50 m	13505
6	3&4	15 m	13601
6	3&4	25 m	13602
6	3&4	50 m	13605

CATGUT

PLAIN

Natural collagen
Twisted multifilament with
a monofilament aspect
Absorbable



Metric EP	USP	Lenght	Order number
2	4/0	100 m	1020
3	3/0	100 m	1030
3,5	2/0	100 m	1035
4	0	100 m	1040
5	1	75 m	1050
6	2	50 m	1060
7	3	50 m	1070
8	4	30 m	1080
9	5	30 m	1090
10	6	30 m	1100

CATGUT

CHROM

Natural collagen
Twisted multifilament with
a monofilament aspect
Absorbable



Metric EP	USP	Lenght	Order number
2	4/0	100 m	2020
3	3/0	100 m	2030
3,5	2/0	100 m	2035
4	0	100 m	2040
5	1	75 m	2050
6	2	50 m	2060
7	3	50 m	2070
8	4	30 m	2080
9	5	30 m	2090
10	6	30 m	2100

SILK

BRAIDED

BLACK

Braided multifilament
Non-absorbable



Metric EP	USP	Lenght	Order number
1,5	4/0	100 m	8015
2	3/0	100 m	8020
3	2/0	100 m	8030
3,5	0	100 m	8035
4	1	100 m	8040
5	2	100 m	8050
6	3 & 4	80 m	8060
7	5	60 m	8070
8	6	50 m	8080
9	7	40 m	8090
10	8	30 m	8100

SUPRAMID

WHITE

Polyamide Twisted and coated
multifilament with a
monofilament aspect
Non-absorbable



Metric EP	USP	Lenght	Order number
1,5	4/0	100 m	3015
2	3/0	100 m	3020
2,5	2/0	100 m	3025
3	2/0	100 m	3030
3,5	0	100 m	3035
4	1	100 m	3040
5	2	75 m	3050
6	3 & 4	50 m	3060
7	5	40 m	3070
8	6	40 m	3080
9	7	30 m	3090

SUPRAMID

BLACK

Polyamide Twisted and
coated multifilament with a
monofilament aspect
Non-absorbable



Metric EP	USP	Lenght	Order number
1,5	4/0	100 m	4015
2	3/0	100 m	4020
2,5	2/0	100 m	4025
3	2/0	100 m	4030
3,5	0	100 m	4035
4	1	100 m	4040
5	2	75 m	4050
6	3 & 4	50 m	4060
7	5	40 m	4070
8	6	40 m	4080
9	7	30 m	4090



LINAFIL

NATURAL

Linen
Twisted multifilament
Non-absorbable



Metric EP	USP	Length	Order number
3	2/0	100 m	5030
3,5	0	100 m	5035
4	1	100 m	5040
5	2	75 m	5050
6	3 & 4	50 m	5060
7	5	40 m	5070

VITAFIL

POLYESTER

WHITE

Polyester
Braided multifilament
Non-absorbable



Metric EP	USP	Length	Order number
2,5	2/0	100 m	6025
3	2/0	100 m	6030
3,5	0	100 m	6035
4	1	100 m	6040
5	2	100 m	6050
6	3 & 4	80 m	6060
7	5	60 m	6070
8	6	50 m	6080
9	7	40 m	6090
10	8	30 m	6100

DERMAFIL

POLYESTER

GREEN

Polyester
Braided multifilament
Non-absorbable



Metric EP	USP	Length	Order number
1,5	4/0	100 m	7015
2	3/0	100 m	7020
2,5	2/0	100 m	7025
3	2/0	100 m	7030
3,5	0	100 m	7035
4	1	100 m	7040
5	2	75 m	7050
6	3 & 4	50 m	7060
7	5	40 m	7070
8	6	40 m	7080
9	7	30 m	7090

The SMI logo is located in the bottom right corner of the page. It consists of the letters 'SMI' in a bold, blue, sans-serif font. The 'S' and 'M' are connected, and the 'I' is separate. The logo is set against a white background with a blue outline.

DACLON

BLUE

Polyamide - Nylon
Monofilament
Non-absorbable



Metric EP	USP	Lenght	Order number
1,5	4/0	100 m	9015
2	3/0	100 m	9020
3	2/0	100 m	9030
3,5	0	100 m	9035
4	1	100 m	9040
5	2	75 m	9050
6	3 & 4	50 m	9060

STEEL

MONOFILAMENT

Non-absorbable



Metric EP	USP	Lenght	Order number
1	5/0	50 m	10010
1,5	4/0	50 m	10015
2	3/0	50 m	10020
3	2/0	50 m	10030
3,5	0	50 m	10035
4	1	50 m	10040