


# DACLON

## NYLON

DACLON  
NYLON



<b>Type</b>	Monofilament
<b>Composition</b>	Extrusion of polyamide 6.0 or 6.6
<b>Coating</b>	None
<b>Colour</b>	Blue Black
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
<b>Presentation</b>	Box with 12 sutures

*For quality  
&  
safety*

**smi**

[www.sutures.be](http://www.sutures.be)

# DACLON

## NYLON

### **CHARACTERISTICS**

- ☑ High tensile strength, non absorbable
- ☑ Synthetic
- ☑ Monofilament with microscopic uniform diameter
- ☑ Pliability and softness
- ☑ Hermetically sealed packing

### **BENEFITS**











- Excellent and permanent tissue support
- Minimal tissue reaction
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
- Good handling properties
- Easy and secure knotting
- Guaranteed seal and product sterility

**INDICATIONS**    General soft tissue closing and/or ligation

**Daclon Blue**  
 especially in    Skin closure  
                       General surgery  
                       Plastic surgery  
                       Gastrointestinal surgery  
                       Gynaecology  
                       Obstetrics  
                       Orthopaedics

**Daclon Black**  
 especially in    Plastic surgery  
                       Ophthalmic surgery

## Round bodied taper point


		Metric	USP	Length	Color	Order number
HR 16 mm		1	5/0	75 cm	black	9100116
		1.5	4/0	75 cm	blue	9150116
HR 20 mm		0.7	6/0	75 cm	black	9070120
		1	5/0	75 cm	black	9100120
		1.5	4/0	75 cm	blue	9150120
		2	3/0	75 cm	blue	9200120
		3	2/0	75 cm	blue	9300120
HR 22 mm		3	2/0	75 cm	blue	9300122
HR 26 mm		1.5	4/0	75 cm	blue	9150126
		2	3/0	75 cm	blue	9200126
		3	2/0	75 cm	blue	9300126
		3.5	0	75 cm	blue	9350126
		4	1	75 cm	blue	9400126
HR 30 mm		3	2/0	75 cm	blue	9300130
		3.5	0	75 cm	blue	9350130
		4	1	75 cm	blue	9400130
HR 37 mm		3	2/0	75 cm	blue	9300137
		3.5	0	75 cm	blue	9350137
HR 40 mm		4	1	75 cm	blue	9400137
		4	1	75 cm	blue	9400140
HR 40 mm		3	2/0	75 cm	blue	9300140
		3.5	0	75 cm	blue	9350140
		4	1	75 cm	blue	9400140
HR 40 mm		5	2	75 cm	blue	9500140
		4	1	75 cm	blue	9400140
		5	2	75 cm	blue	9500140
HR 50 mm		4	1	100 cm	blue	9400150

1/2 Circle




## Round bodied taper point loop

		Metric	USP	Length	Order number
--	--	--------	-----	--------	--------------

HR 30 mm 	3.5	0	150 cm	<i>blue</i>	9350930
	4	1	150 cm	<i>blue</i>	9400930

HR 35 mm 	3.5	0	150 cm	<i>blue</i>	9350935
	4	1	150 cm	<i>blue</i>	9400935

## 1/2 Circle

HR 40 mm 	3.5	0	150 cm	<i>blue</i>	9350940
	4	1	150 cm	<i>blue</i>	9400940

HR 50 mm 	3.5	0	150 cm	<i>blue</i>	9350950
	4	1	150 cm	<i>blue</i>	9400950

## Round bodied taper point double armed

	Metric	USP	Length	Order number
--	--------	-----	--------	--------------


## 1/2 Circle

2 x HR 65 mm 	3	2/0	90 cm	<i>blue</i>	9300165-2
--	---	-----	-------	-------------	-----------

Reverse cutting



3/8 Circle













		Metric	USP	Length		Order number
DS 12 mm		0.5	7/0	75 cm	<i>black</i>	9051512
		0.7	6/0	75 cm	<i>black</i>	9071512
		1	5/0	75 cm	<i>black</i>	9101512
		1.5	4/0	75 cm	<i>blue</i>	9151512
DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>black</i>	9071512P
		1	5/0	50 cm	<i>black</i>	9101512P
		1.5	4/0	50 cm	<i>blue</i>	9151512P
DS 16 mm		0.7	6/0	45 cm	<i>black</i>	9071516-45
		0.7	6/0	75 cm	<i>black</i>	9071516
		1	5/0	45 cm	<i>black</i>	9101516-45
		1	5/0	75 cm	<i>black</i>	9101516
		1.5	4/0	45 cm	<i>blue</i>	9151516-45
		1.5	4/0	75 cm	<i>blue</i>	9151516
DS 16 mm Easy Pass		0.7	6/0	50 cm	<i>black</i>	9071516P
		1	5/0	50 cm	<i>black</i>	9101516P
		1.5	4/0	50 cm	<i>blue</i>	9151516P
DS 16 mm black		1	5/0	75 cm	<i>black</i>	9101516VB
DS 19 mm		1	5/0	45 cm	<i>black</i>	9101519-45
		1	5/0	75 cm	<i>black</i>	9101519
		1.5	4/0	45 cm	<i>blue</i>	9151519-45
		1.5	4/0	75 cm	<i>blue</i>	9151519
		2	3/0	45 cm	<i>blue</i>	9201519-45
		2	3/0	75 cm	<i>blue</i>	9201519
DS 19 mm Easy Pass		1.5	4/0	50 cm	<i>blue</i>	9151519P
		2	3/0	50 cm	<i>blue</i>	9201519P
		3	2/0	50 cm	<i>blue</i>	9301519P
DS 19 mm black		1	5/0	75 cm	<i>black</i>	9101519VB
		1.5	4/0	75 cm	<i>blue</i>	9151519VB

Reverse cutting



3/8 Circle



		Metric	USP	Length		Order number
DS 24 mm		1	5/0	75 cm	<i>black</i>	9101524
		1.5	4/0	45 cm	<i>blue</i>	9151524-45
		1.5	4/0	75 cm	<i>blue</i>	9151524
		2	3/0	45 cm	<i>blue</i>	9201524-45
		2	3/0	75 cm	<i>blue</i>	9201524
		3	2/0	45 cm	<i>blue</i>	9301524-45
		3	2/0	75 cm	<i>blue</i>	9301524
		3.5	0	75 cm	<i>blue</i>	9351524
		4	1	75 cm	<i>blue</i>	9401524
		5	2	90 cm	<i>blue</i>	9501524
DS 24 mm Easy Pass		3	2/0	50 cm	<i>blue</i>	9301524P
DS 26 mm		2	3/0	45 cm	<i>blue</i>	9201526-45
		2	3/0	75 cm	<i>blue</i>	9201526
		3	2/0	45 cm	<i>blue</i>	9301526-45
		3	2/0	75 cm	<i>blue</i>	9301526
DS 30 mm		2	3/0	75 cm	<i>blue</i>	9201530
		3	2/0	75 cm	<i>blue</i>	9301530
		3.5	0	75 cm	<i>blue</i>	9351530
		4	1	75 cm	<i>blue</i>	9401530
DS 36 mm		2	3/0	75 cm	<i>blue</i>	9201536
		3	2/0	75 cm	<i>blue</i>	9301536
		3.5	0	75 cm	<i>blue</i>	9351536
DS 40 mm		3	2/0	75 cm	<i>blue</i>	9301540
		3.5	0	75 cm	<i>blue</i>	9351540
		4	1	75 cm	<i>blue</i>	9401540
DS 45 mm		5	2	90 cm	<i>blue</i>	9501540
		3	2/0	100 cm	<i>blue</i>	9301545
		4	1	75 cm	<i>blue</i>	9401545
DS 60 mm		3	2/0	75 cm	<i>blue</i>	9301560
		3.5	0	75 cm	<i>blue</i>	9351560
DS 90 mm		3	2/0	100 cm	<i>blue</i>	9301590
		3.5	0	100 cm	<i>blue</i>	9351590
		4	1	100 cm	<i>blue</i>	9401590

## Reverse cutting



		Metric	USP	Length	Order number
1/2 Circle 	HS 17 mm	1.5	4/0	75 cm	<i>blue</i> 9151417
		2	3/0	75 cm	<i>blue</i> 9201417
	HS 22 mm	2	3/0	75 cm	<i>blue</i> 9201422
		HS 26 mm	2	3/0	75 cm
	3		2/0	75 cm	<i>blue</i> 9301426
	HS 30 mm	3.5	0	75 cm	<i>blue</i> 9351430
		HS 36 mm	3.5	0	75 cm
	4		1	75 cm	<i>blue</i> 9401436
	5		2	75 cm	<i>blue</i> 9501436
	HS 40 mm	3.5	0	75 cm	<i>blue</i> 9351440
4		1	75 cm	<i>blue</i> 9401440	
HS 48 mm	4	1	100 cm	<i>blue</i> 9401448	

## Cutting
















		Metric	USP	Length	Order number
Straight	GS 60 mm	2	3/0	75 cm	<i>blue</i> 9202160
		3	2/0	75 cm	<i>blue</i> 9302160






## Cutting



		Metric	USP	Length	Order number
3/8 Circle 	NS 26 mm	2	3/0	75 cm	<i>blue</i> 9201726
	NS 30 mm	3	2/0	75 cm	<i>blue</i> 9301730
	NS 40 mm	3	2/0	75 cm	<i>blue</i> 9301740

Micro needle taper point			Metric	USP	Length	Order number	
3/8 Circle 	MDR 3,3 mm 50µm		0.1	11/0	10 cm	<i>black</i>	MN111
	MDR 3,8 mm 50µm		0.2	10/0	10 cm	<i>black</i>	MN102
	MDR 5,1 mm 70µm		0.01	12/0	10 cm	<i>black</i>	MN124
			0.1	11/0	10 cm	<i>black</i>	MN114
			0.2	10/0	10 cm	<i>black</i>	MN104
	MDR 3,8 mm 70µm		0.1	11/0	10 cm	<i>black</i>	MN113
	MDR 3,8 mm 100µm		0.2	10/0	15 cm	<i>black</i>	MN105
	MDR 5,1 mm 100µm		0.2	10/0	15 cm	<i>black</i>	MN106
0.3			9/0	15 cm	<i>black</i>	MN096	
MDR 6,4 mm 140µm		0.4	8/0	15 cm	<i>black</i>	MN087	
MDR 9,7 mm 200µm		0.5	7/0	45 cm	<i>black</i>	MN078	

Round bodied taper point double armed			Metric	USP	Length	Order number	
1/4 Circle 	2 x VR 8 mm 220 µm		0.4	8/0	45 cm	<i>black</i>	ON81

Side cutting spatula double armed			Metric	USP	Length	Order number	
1/2 Circle 	2 x HSP 6,2 mm 150 µm		0.2	10/0	30 cm	<i>black</i>	ON102
	2 x HSP 6,4 mm 220 µm		0.2	10/0	30 cm	<i>black</i>	ON106
	2 x HSP 8 mm 300 µm		0.7	6/0	45 cm	<i>black</i>	ON62



## Side cutting spatula double armed



### 3/8 Circle




			Metric	USP	Length		Order number
3/8 Circle	2 x DSP 6 mm 220 µm		0.2	10/0	30 cm	<i>black</i>	ON103
			0.1	11/0	20 cm	<i>black</i>	ON111
	2 x DSP 6,2 mm 150 µm		0.2	10/0	30 cm	<i>black</i>	ON101
			0.3	9/0	30 cm	<i>black</i>	ON95
			0.3	9/0	45 cm	<i>black</i>	ON95-45
			0.4	8/0	30 cm	<i>black</i>	ON85
			0.4	8/0	45 cm	<i>black</i>	ON85-45
			0.5	7/0	30 cm	<i>black</i>	ON75
	2 x DSP 6,4 mm 220 µm		0.2	10/0	30 cm	<i>black</i>	ON105
			0.2	10/0	45 cm	<i>black</i>	ON105-45
2 x DSP 8 mm 300 µm		0.7	6/0	45 cm	<i>black</i>	ON61	

## Side cutting spatula double armed



### 1/4 Circle



			Metric	USP	Length		Order number
1/4 Circle	2 x VSP 8 mm 350 µm		1	5/0	45 cm	<i>black</i>	ON51