


# DACLON NYLON

DACLON  
NYLON



<b>Type</b>	Monofilament
<b>Composition</b>	Extrusion of polyamide 6.0 or 6.6
<b>Coating</b>	None
<b>Colour</b>	Blue Black
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
<b>Presentation</b>	Box with 12 sutures

*For quality  
&  
safety*

**smi**

[www.sutures.be](http://www.sutures.be)

# DACLON

## NYLON

### **CHARACTERISTICS**

- ☑ High tensile strength, non absorbable
- ☑ Synthetic
- ☑ Monofilament with microscopic uniform diameter
- ☑ Pliability and softness
- ☑ Hermetically sealed packing

### **BENEFITS**

- Excellent and permanent tissue support
- Minimal tissue reaction
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
- Good handling properties
- Easy and secure knotting
- Guaranteed seal and product sterility




**INDICATIONS**    General soft tissue closing and/or ligation

**Daclon Blue**  
 especially in    Skin closure  
                       General surgery  
                       Plastic surgery  
                       Gastrointestinal surgery  
                       Gynaecology  
                       Obstetrics  
                       Orthopaedics

**Daclon Black**  
 especially in    Plastic surgery  
                       Ophthalmic surgery







## Round bodied taper point







### 1/2 Circle






		Metric	USP	Length	Color	Order number
HR 16 mm		1	5/0	75 cm	black	9100116
		1.5	4/0	75 cm	blue	9150116
HR 20 mm		0.7	6/0	75 cm	black	9070120
		1	5/0	75 cm	black	9100120
		1.5	4/0	75 cm	blue	9150120
		2	3/0	75 cm	blue	9200120
		3	2/0	75 cm	blue	9300120
HR 22 mm		3	2/0	75 cm	blue	9300122

## Reverse cutting

### 3/8 Circle

		Metric	USP	Length	Color	Order number
DS 12 mm		0.5	7/0	75 cm	black	9051512
		0.7	6/0	75 cm	black	9071512
		1	5/0	75 cm	black	9101512
		1.5	4/0	75 cm	blue	9151512
DS 12 mm Easy Pass		0.7	6/0	50 cm	black	9071512P
		1	5/0	50 cm	black	9101512P
		1.5	4/0	50 cm	blue	9151512P
DS 16 mm		0.7	6/0	45 cm	black	9071516-45
		0.7	6/0	75 cm	black	9071516
		1	5/0	45 cm	black	9101516-45
		1	5/0	75 cm	black	9101516
		1.5	4/0	45 cm	blue	9151516-45
		1.5	4/0	75 cm	blue	9151516
		2	3/0	75 cm	blue	9201516
DS 16 mm Easy Pass		0.7	6/0	50 cm	black	9071516P
		1	5/0	50 cm	black	9101516P
		1.5	4/0	50 cm	blue	9151516P
DS 16 mm black		1	5/0	75 cm	black	9101516VB
DS 19 mm		1	5/0	45 cm	black	9101519-45
		1	5/0	75 cm	black	9101519
		1.5	4/0	45 cm	blue	9151519-45
		1.5	4/0	75 cm	blue	9151519
		2	3/0	45 cm	blue	9201519-45
		2	3/0	75 cm	blue	9201519
		3	2/0	75 cm	blue	9301519

Reverse cutting 		Metric	USP	Length	Order number
<b>3/8 Circle</b> 	DS 19 mm Easy Pass 	1.5	4/0	50 cm	<i>blue</i> 9151519P
		2	3/0	50 cm	<i>blue</i> 9201519P
		3	2/0	50 cm	<i>blue</i> 9301519P
	DS 19 mm black 	1	5/0	75 cm	<i>black</i> 9101519VB
		1.5	4/0	75 cm	<i>blue</i> 9151519VB
	DS 24 mm 	1	5/0	75 cm	<i>black</i> 9101524
		1.5	4/0	45 cm	<i>blue</i> 9151524-45
		1.5	4/0	75 cm	<i>blue</i> 9151524
		2	3/0	45 cm	<i>blue</i> 9201524-45
		2	3/0	75 cm	<i>blue</i> 9201524
		3	2/0	45 cm	<i>blue</i> 9301524-45
	DS 24 mm Easy Pass 	3	2/0	75 cm	<i>blue</i> 9301524
DS 24 mm Easy Pass	3	2/0	50 cm	<i>blue</i> 9301524P	

Reverse cutting 		Metric	USP	Length	Order number
<b>1/2 Circle</b> 	HS 17 mm 	1.5	4/0	75 cm	<i>blue</i> 9151417
		2	3/0	75 cm	<i>blue</i> 9201417
	HS 22 mm 	2	3/0	75 cm	<i>blue</i> 9201422
	HS 26 mm 	2	3/0	75 cm	<i>blue</i> 9201426
		3	2/0	75 cm	<i>blue</i> 9301426