

POLYESTER



Type	Braided multifilament
Composition	Polyester - a polymer of polyethylene terephthalate
Coating	Silicone
Colour	Green
Tissue reaction	Minimal
Absorption	Non absorbable, gradually encapsulated by connective tissue
Presentation	Box with 12 sutures

*For quality
&
safety*

smi

www.sutures.be

POLYESTER

CHARACTERISTICS

- ✓ Unsurpassed tensile strength, non absorbable
- ✓ Resistance to repeated bending
- ✓ Synthetic
- ✓ Braided structure

- ✓ Coated multifilament

- ✓ Hermetically sealed packing

INDICATIONS General soft tissue closing and/or ligation
especially in Cardiovascular surgery
Ophthalmic surgery
General surgery
Skin closure
Orthopaedics
Gastrointestinal surgery
Obstetrics

BENEFITS

- Excellent and permanent tissue support






- Use in cardiovascular surgery (heart prosthesis)




- Minimal tissue reaction




- Excellent handling properties
High flexibility
High tensile strength
Excellent knot security







- Soft passage through the tissues with minimal sawing, tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action

- Guaranteed seal and product sterility







 Round bodied taper point	Metric	USP	Length	Order number			
1/2 circle 	HR 16 mm		1	5/0	75 cm	<i>green</i>	7100116
		1.5	4/0	75 cm	<i>green</i>	7150116	
		2	3/0	75 cm	<i>green</i>	7200116	
	HR 20 mm		1.5	4/0	75 cm	<i>green</i>	7150120
		2	3/0	75 cm	<i>green</i>	7200120	
		2	3/0	75 cm	<i>green</i>	7200126	
HR 26 mm		3	2/0	75 cm	<i>green</i>	7300126	

 Round bodied taper point	Metric	USP	Length	Order number			
Straight 	GR 19 mm		1.5	4/0	45 cm	<i>green</i>	7152019
		2	3/0	50 cm	<i>green</i>	7202019	

 Round bodied taper cutting point	Metric	USP	Length	Order number			
1/2 circle 	HRT 18 mm		1.5	4/0	75 cm	<i>green</i>	7150518
		2	3/0	75 cm	<i>green</i>	7200518	

 Reverse cutting	Metric	USP	Length	Order number			
3/8 circle 	DS 12 mm		1	5/0	45 cm	<i>green</i>	7101512
		1	5/0	75 cm	<i>green</i>	7101516	
	DS 16 mm		1.5	4/0	75 cm	<i>green</i>	7151516
		2	3/0	75 cm	<i>green</i>	7201516	
	DS 19 mm		1	5/0	75 cm	<i>green</i>	7101519
		1.5	4/0	75 cm	<i>green</i>	7151519	
		2	3/0	75 cm	<i>green</i>	7201519	
	DS 19 mm black		3	2/0	75 cm	<i>green</i>	7301519
		1.5	4/0	75 cm	<i>green</i>	7151519VB	
	DS 24 mm	2	3/0	75 cm	<i>green</i>	7201519VB	
		2	3/0	75 cm	<i>green</i>	7201524	
		3	2/0	75 cm	<i>green</i>	7301524	

POLYESTER

 Reverse cutting		Metric	USP	Length		Order number	
	HS 17 mm		1.5	4/0	75 cm	<i>green</i>	7151417
			2	3/0	75 cm	<i>green</i>	7201417
1/2 circle	HS 20 mm		3	2/0	75 cm	<i>green</i>	7301420
	HS 22 mm		2	3/0	75 cm	<i>green</i>	7201422
	HS 26 mm		2	3/0	75 cm	<i>green</i>	7201426
			3	2/0	75 cm	<i>green</i>	7301426