


POLYPROPYLENE

BLUE



Type	Monofilament
Composition	Polypropylene, a polymer of propylene
Coating	None
Colour	Blue
Tissue reaction	Minimal
Absorption	Non absorbable
Presentation	Box with 12 sutures

*For quality
&
safety*



smi

www.sutures.be

POLYPROPYLENE

BLUE







CHARACTERISTICS




- ✓ High tensile strength, non absorbable
- ✓ Resistance to repeated bending
- ✓ Synthetic
- ✓ Monofilament with microscopic uniform diameter
- ✓ Pliability and softness
- ✓ Hermetically sealed packing









INDICATIONS General soft tissue closing and/or ligation
especially in Cardiovascular surgery
Neurosurgery
Ophthalmic surgery
Microsurgery
Plastic surgery
Skin closure
Orthopaedics
Hernias (also as a mesh)
Gastrointestinal surgery
General surgery
Gynaecology
Obstetrics








BENEFITS





Excellent and permanent tissue support
Use in cardiovascular surgery (heart prosthesis or vascular anastomosis)
Minimal tissue reaction
Soft passage through the tissues without sawing, tissue drag and trauma
Absence of capillarity (no support of bacterial growth)
Excellent knot tie down
Good handling properties
Easy and secure knotting
Guaranteed seal and product sterility




Round bodied taper point 		Metric	USP	Length	Order number
1/2 Circle 	HR 10 mm 	0.5	7/0	45 cm	<i>blue</i> 5050110
		0.7	6/0	75 cm	<i>blue</i> 5070110
	HR 15 mm 	0.7	6/0	75 cm	<i>blue</i> 5070115
		1	5/0	75 cm	<i>blue</i> 5100115
	HR 17 mm 	1.5	4/0	75 cm	<i>blue</i> 5150115
		1	5/0	75 cm	<i>blue</i> 5100117
	HR 20 mm 	2	3/0	75 cm	<i>blue</i> 5200117
		1	5/0	75 cm	<i>blue</i> 5100120
		1.5	4/0	75 cm	<i>blue</i> 5150120
		2	3/0	75 cm	<i>blue</i> 5200120

Round bodied taper point 		Metric	USP	Length	Order number
3/8 Circle 	DR 16 mm 	1.5	4/0	75 cm	<i>blue</i> 5150216

Reverse cutting 		Metric	USP	Length	Order number
3/8 Circle 	DS 10 mm 	0.7	6/0	45 cm	<i>blue</i> 5071510
		0.5	7/0	45 cm	<i>blue</i> 5051512
	DS 12 mm 	0.7	6/0	45 cm	<i>blue</i> 5071512
		1	5/0	45 cm	<i>blue</i> 5101512
		1.5	4/0	75 cm	<i>blue</i> 5151512
	DS 12 mm Easy Pass 	0.7	6/0	50 cm	<i>blue</i> 5071512P
		1	5/0	50 cm	<i>blue</i> 5101512P
	DS 16 mm 	0.7	6/0	45 cm	<i>blue</i> 5071516
		1	5/0	45 cm	<i>blue</i> 5101516
		1.5	4/0	75 cm	<i>blue</i> 5151516
		2	3/0	75 cm	<i>blue</i> 5201516
		0.7	6/0	50 cm	<i>blue</i> 5071516P
	DS 16 mm Easy Pass 	1	5/0	50 cm	<i>blue</i> 5101516P
		1.5	4/0	50 cm	<i>blue</i> 5151516P
	DS 16 mm black 	1	5/0	45 cm	<i>blue</i> 5101516VB

Reverse cutting		Metric	USP	Length	Order number	
3/8 Circle 		0.7	6/0	75 cm	<i>blue</i> 5071519	
		1	5/0	75 cm	<i>blue</i> 5101519	
		1.5	4/0	45 cm	<i>blue</i> 5151519-45	
		1.5	4/0	75 cm	<i>blue</i> 5151519	
		2	3/0	45 cm	<i>blue</i> 5201519-45	
		2	3/0	75 cm	<i>blue</i> 5201519	
		3	2/0	75 cm	<i>blue</i> 5301519	
		DS 19 mm Easy Pass		1	5/0	50 cm
	1.5			4/0	50 cm	<i>blue</i> 5151519P
	2			3/0	50 cm	<i>blue</i> 5201519P
	DS 19 mm black		1	5/0	75 cm	<i>blue</i> 5101519VB
			1.5	4/0	75 cm	<i>blue</i> 5151519VB
	DS 24 mm		1.5	4/0	75 cm	<i>blue</i> 5151524
			2	3/0	45 cm	<i>blue</i> 5201524-45
			2	3/0	75 cm	<i>blue</i> 5201524
			3	2/0	45 cm	<i>blue</i> 5301524-45
			3	2/0	75 cm	<i>blue</i> 5301524
	DS 24 mm Easy Pass		2	3/0	50 cm	<i>blue</i> 5201524P
3			2/0	50 cm	<i>blue</i> 5301524P	

Reverse cutting		Metric	USP	Length	Order number	
1/2 Circle 		HS 22 mm	2	3/0	45 cm	<i>blue</i> 5201422
		HS 24 mm	3	2/0	75 cm	<i>blue</i> 5301424

Round bodied taper cutting point		Metric	USP	Length	Order number	
1/2 Circle 		HRT 10 mm	0.7	6/0	75 cm	<i>blue</i> 5070510