


SILK



	<i>Silk Braided</i>	<i>Silk Virgin</i>
Type	Braided multifilament	Twisted multifilament
Composition	Braided fibres from the cocoon of the silkworm	Twisted fibres from the cocoon of the silkworm
Coating	Waxes Silicone	None
Colour	Black	Blue
Tissue reaction	Moderate	
Absorption	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.	
Presentation	Box with 12 sutures	

*For quality
&
safety*

smi

www.sutures.be

SILK

CHARACTERISTICS

- ☑ High tensile strength, non absorbable
- ☑ Braided or twisted structure

- ☑ Coated multifilament

- ☑ Hermetically sealed packing

BENEFITS

- Good and extended tissue support for up to 3 month

- Excellent handling properties
High flexibility
High tensile strength
Excellent knot security

- Soft passage through the tissues with minimal sawing, tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation
especially in General surgery
Skin closure
Ophthalmic surgery
Oral surgery
Gastrointestinal surgery
Gynaecology
Obstetrics

CONTRAINDICATIONS Urinary tract
Biliary tract

Known allergies or sensitivities to silk

Round bodied taper point



1/2 Circle



		Metric	USP	Length	Color	Order number	
HR 12 mm		0.7	6/0	45 cm	black	8070112	
		1	5/0	45 cm	black	8100112	
		1.5	4/0	75 cm	black	8150112	
	HR 16 mm		0.7	6/0	45 cm	black	8070116
			1	5/0	75 cm	black	8100116
			1.5	4/0	75 cm	black	8150116
HR 20 mm		2	3/0	75 cm	black	8200116	
		1	5/0	75 cm	black	8100120	
		1.5	4/0	75 cm	black	8150120	
		2	3/0	75 cm	black	8200120	
		3	2/0	75 cm	black	8300120	

Round bodied taper point



3/8 Circle



		Metric	USP	Length	Color	Order number
DR 12 mm		1	5/0	45 cm	black	8100212
		1.5	4/0	75 cm	black	8150216
		2	3/0	75 cm	black	8200216

Reverse cutting



3/8 Circle



		Metric	USP	Length	Color	Order number
DS 12 mm		0.5	7/0	45 cm	black	8051512
		0.7	6/0	45 cm	black	8071512
		1	5/0	45 cm	black	8101512
		1.5	4/0	75 cm	black	8151512
		2	3/0	75 cm	black	8201512
DS 12 mm Easy Pass		0.7	6/0	50 cm	black	8071512P
DS 14 mm		0.7	6/0	45 cm	black	8071514
DS 16 mm		0.7	6/0	75 cm	black	8071516
		1	5/0	75 cm	black	8101516
		1.5	4/0	75 cm	black	8151516
		2	3/0	75 cm	black	8201516
		3	2/0	75 cm	black	8301516
DS 16 mm Easy Pass		1	5/0	50 cm	black	8101516P
		1.5	4/0	50 cm	black	8151516P

Reverse cutting










		Metric	USP	Length	Order number
3/8 Circle 	DS 16 mm black	1	5/0	75 cm	<i>black</i> 8101516VB
		1.5	4/0	75 cm	<i>black</i> 8151516VB
	DS 19 mm	0.7	6/0	75 cm	<i>black</i> 8071519
		1	5/0	75 cm	<i>black</i> 8101519
		1.5	4/0	75 cm	<i>black</i> 8151519
		2	3/0	45 cm	<i>black</i> 8201519-45
		2	3/0	75 cm	<i>black</i> 8201519
		3	2/0	75 cm	<i>black</i> 8301519
	DS 19 mm Easy Pass	1.5	4/0	50 cm	<i>black</i> 8151519P
		2	3/0	50 cm	<i>black</i> 8201519P
	DS 19 mm black	1	5/0	75 cm	<i>black</i> 8101519VB
		1.5	4/0	75 cm	<i>black</i> 8151519VB
	DS 24 mm	2	3/0	75 cm	<i>black</i> 8201519VB
		1	5/0	75 cm	<i>black</i> 8101524
		1.5	4/0	75 cm	<i>black</i> 8151524
	DS 24 mm Easy Pass	2	3/0	75 cm	<i>black</i> 8201524
		3	2/0	75 cm	<i>black</i> 8301524
		2	3/0	50 cm	<i>black</i> 8201524P

Reverse cutting



		Metric	USP	Length	Order number
1/2 Circle 	HS 12 mm	0.7	6/0	75 cm	<i>black</i> 8071412
	HS 17 mm	1	5/0	75 cm	<i>black</i> 8101417
		1.5	4/0	75 cm	<i>black</i> 8151417
		2	3/0	75 cm	<i>black</i> 8201417
	HS 20 mm	1.5	4/0	75 cm	<i>black</i> 8151420
		2	3/0	75 cm	<i>black</i> 8201420
	HS 22 mm	1	5/0	75 cm	<i>black</i> 8101422
		1.5	4/0	75 cm	<i>black</i> 8151422
		2	3/0	75 cm	<i>black</i> 8201422
	HS 26 mm	3	2/0	75 cm	<i>black</i> 8301422
		1	5/0	75 cm	<i>black</i> 8101426
		1.5	4/0	75 cm	<i>black</i> 8151426
	HS 26 mm	2	3/0	75 cm	<i>black</i> 8201426
		3	2/0	75 cm	<i>black</i> 8301426

Round bodied taper cutting point 		Metric	USP	Length	Order number		
1/2 Circle 	HRT 18 mm		1	5/0	75 cm	<i>black</i>	8100518
			1.5	4/0	75 cm	<i>black</i>	8150518
			2	3/0	75 cm	<i>black</i>	8200518
	HRT 27 mm		3	2/0	75 cm	<i>black</i>	8300518
			2	3/0	75 cm	<i>black</i>	8200527
			3	2/0	75 cm	<i>black</i>	8300527

Round bodied taper cutting point 		Metric	USP	Length	Order number		
3/8 Circle 	DRT 18 mm		1.5	4/0	75 cm	<i>black</i>	8150418
			2	3/0	75 cm	<i>black</i>	8200418