


SUPRAMID

BLACK



Type	Twisted and coated multifilament Sizes 5/0 and 6/0 monofilament
Composition	Polyamide 6.6 strands enclosed in a polyamide 6 sheet
Colour	Black
Tissue reaction	Minimal
Absorption	Non absorbable, gradually encapsulated by connective tissue
Presentation	Box with 12 sutures

*For quality
&
safety*

smi

www.sutures.be

SUPRAMID

BLACK

CHARACTERISTICS

- ✓ High tensile strength, non absorbable
- ✓ Synthetic
- ✓ Twisted structure

- ✓ Multifilament in a sheath

- ✓ Pliability and softness

- ✓ Hermetically sealed packing

INDICATIONS General soft tissue closing
 and/or ligation
especially in Skin closure

BENEFITS

Excellent and permanent tissue support





Minimal tissue reaction




Excellent handling properties
High flexibility
High tensile strength
Excellent knot security







Soft passage through the tissues without sawing,
tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action




Good handling properties
Easy and secure knotting










Guaranteed seal and product sterility




Round bodied taper point			Metric	USP	Length	Order number
1/2 Circle 	HR 16 mm		1.5	4/0	75 cm	<i>black</i> 4150116
			2	3/0	75 cm	<i>black</i> 4200116
	HR 20 mm		2	3/0	45 cm	<i>black</i> 4200120




Round bodied taper point			Metric	USP	Length	Order number
Straight 	GR 22 mm		1.5	4/0	50 cm	<i>black</i> 4152022

Reverse cutting			Metric	USP	Length	Order number
1/2 Circle 	HS 17 mm		1.5	4/0	75 cm	<i>black</i> 4151417
			2	3/0	75 cm	<i>black</i> 4201417
	HS 20 mm		1	5/0	75 cm	<i>black</i> 4101420
			1.5	4/0	75 cm	<i>black</i> 4151420
			2	3/0	75 cm	<i>black</i> 4201420
	HS 22 mm		1.5	4/0	75 cm	<i>black</i> 4151422
HS 26 mm		1.5	4/0	75 cm	<i>black</i> 4151426	
		2	3/0	75 cm	<i>black</i> 4201426	
		3	2/0	75 cm	<i>black</i> 4301426	

Cutting			Metric	USP	Length	Order number
Straight 	GS 20 mm		1.5	4/0	50 cm	<i>black</i> 4152120

Reverse cutting			Metric	USP	Length	Order number
3/8 Circle 	DS 12 mm		0.7	6/0	75 cm	<i>black</i> 4071512
			1	5/0	75 cm	<i>black</i> 4101512
			1.5	4/0	45 cm	<i>black</i> 4151512
	DS 16 mm		0.7	6/0	75 cm	<i>black</i> 4071516
			1	5/0	75 cm	<i>black</i> 4101516
			1.5	4/0	45 cm	<i>black</i> 4151516
			1.5	4/0	75 cm	<i>black</i> 4151516-75
			2	3/0	45 cm	<i>black</i> 4201516
			2	3/0	75 cm	<i>black</i> 4201516-75
	DS 16 mm black		1	5/0	75 cm	<i>black</i> 4101516VB
	DS 19 mm		0.7	6/0	75 cm	<i>black</i> 4071519
			1	5/0	75 cm	<i>black</i> 4101519
			1.5	4/0	45 cm	<i>black</i> 4151519
			1.5	4/0	75 cm	<i>black</i> 4151519-75
			2	3/0	75 cm	<i>black</i> 4201519
			3	2/0	75 cm	<i>black</i> 4301519
	DS 19 mm Easy Pass		1.5	4/0	75 cm	<i>black</i> 4151519P-75
	DS 19 mm black		1	5/0	75 cm	<i>black</i> 4101519VB
1.5			4/0	75 cm	<i>black</i> 4151519VB	
DS 24 mm		2	3/0	75 cm	<i>black</i> 4201524	
		3	2/0	75 cm	<i>black</i> 4301524	

Round bodied taper cutting point			Metric	USP	Length	Order number
1/2 Circle 	HRT 18 mm		1.5	4/0	75 cm	<i>black</i> 4150518
			2	3/0	75 cm	<i>black</i> 4200518

Round bodied taper cutting point			Metric	USP	Length	Order number
3/8 Circle 	DRT 18 mm		1.5	4/0	75 cm	<i>black</i> 4150418
			2	3/0	75 cm	<i>black</i> 4200418