

SURGICRYL[®]


PGA

POLYGLYCOLIC ACID

SURGICRYL[®]
PGA
POLYGLYCOLIC ACID



Type	Braided and coated multifilament
Composition	Polyglycolic acid
Coating	Polycaprolactone and calcium stearate (1%)
Colour	Violet / Beige - undyed
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days. Approximately 50% of tensile strength remain after 21 days.
Presentation	Box with 12 sutures



*For quality
&
safety*

smi

www.sutures.be

SURGICRYL[®]

PGA
POLYGLYCOLIC ACID

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure

- Coated multifilament

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
 - High flexibility
 - High tensile strength
 - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Gynaecology
- Obstetrics
- Plastic surgery
- Urology
- Ophthalmic surgery
- Orthopaedics

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an extended wound support is required.






Caution in patients with delayed wound healing (e.g. compromised immune system).




TENSILE STRENGTH










ABSORPTION



Round bodied taper point 		Metric	USP	Length	Order number	
1/2 Circle 	HR 13 mm 	0.7	6/0	75 cm	<i>violet</i>	11070113
		1	5/0	75 cm	<i>violet</i>	11100113
		1.5	4/0	75 cm	<i>violet</i>	11150113
	HR 17 mm 	0.7	6/0	75 cm	<i>violet</i>	11070117
		1	5/0	75 cm	<i>violet</i>	11100117
		1.5	4/0	75 cm	<i>violet</i>	11150117
		2	3/0	75 cm	<i>violet</i>	11200117
		3	2/0	75 cm	<i>violet</i>	11300117
		HR 22 mm 	1.5	4/0	75 cm	<i>violet</i>
	2		3/0	75 cm	<i>violet</i>	11200122
	2		3/0	75 cm	<i>undyed</i>	11200122U
	3		2/0	75 cm	<i>violet</i>	11300122

Round bodied taper point 		Metric	USP	Length	Order number	
3/8 Circle 	DR 17 mm 	1.5	4/0	75 cm	<i>violet</i>	11150217
		2	3/0	75 cm	<i>violet</i>	11200217

Reverse cutting 		Metric	USP	Length	Order number	
1/2 Circle 	HS 22 mm 	2	3/0	75 cm	<i>violet</i>	11201422
		2	3/0	75 cm	<i>undyed</i>	11201422U
		3	2/0	75 cm	<i>violet</i>	11301422
		3	2/0	75 cm	<i>undyed</i>	11301422U
	HS 26 mm 	1.5	4/0	75 cm	<i>violet</i>	11151426
		2	3/0	75 cm	<i>violet</i>	11201426
		3	2/0	75 cm	<i>violet</i>	11301426

Cutting 		Metric	USP	Length	Order number	
3/8 Circle 	NS 16 mm 	1	5/0	45 cm	<i>violet</i>	11101716

Reverse cutting



3/8 Circle



		Metric	USP	Length		Order number
DS 12 mm		0.7	6/0	45 cm	<i>violet</i>	11071512
		0.7	6/0	45 cm	<i>undyed</i>	11071512U
		1	5/0	45 cm	<i>violet</i>	11101512
		1.5	4/0	75 cm	<i>violet</i>	11151512
DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>undyed</i>	11071512UP
		1	5/0	50 cm	<i>undyed</i>	11101512UP
		1.5	4/0	50 cm	<i>undyed</i>	11151512UP
DS 16 mm		0.7	6/0	75 cm	<i>violet</i>	11071516
		1	5/0	75 cm	<i>violet</i>	11101516
		1	5/0	75 cm	<i>undyed</i>	11101516U
		1.5	4/0	75 cm	<i>violet</i>	11151516
		1.5	4/0	75 cm	<i>undyed</i>	11151516U
		2	3/0	75 cm	<i>violet</i>	11201516
DS 16 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i>	11151516UP
DS 19 mm		1	5/0	75 cm	<i>violet</i>	11101519
		1.5	4/0	75 cm	<i>violet</i>	11151519
		1.5	4/0	75 cm	<i>undyed</i>	11151519U
		2	3/0	75 cm	<i>violet</i>	11201519
		2	3/0	75 cm	<i>undyed</i>	11201519U
		3	2/0	75 cm	<i>violet</i>	11301519
DS 19 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i>	11151519UP
		2	3/0	50 cm	<i>undyed</i>	11201519UP
DS 24 mm		2	3/0	75 cm	<i>violet</i>	11201524
		2	3/0	75 cm	<i>undyed</i>	11201524U
		3	2/0	75 cm	<i>violet</i>	11301524
		3	2/0	75 cm	<i>undyed</i>	11301524U
DS 24 mm Easy Pass		2	3/0	50 cm	<i>undyed</i>	11201524UP
		3	2/0	50 cm	<i>undyed</i>	11301524UP