


# DACLON

## NYLON

DACLON  
NYLON



<b>Type</b>	Monofilament
<b>Composition</b>	Extrusion of polyamide 6.0 or 6.6
<b>Coating</b>	None
<b>Colour</b>	Black
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
<b>Presentation</b>	Box with 12 sutures

*For quality  
&  
safety*

**smi**

[www.sutures.be](http://www.sutures.be)

# DACLON

## NYLON

### **CHARACTERISTICS**




- ☑ High tensile strength, non absorbable
- ☑ Synthetic
- ☑ Monofilament with microscopic uniform diameter
- ☑ Pliability and softness
- ☑ Hermetically sealed packing




### **BENEFITS**







- Excellent and permanent tissue support
- Minimal tissue reaction
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
- Good handling properties
- Easy and secure knotting
- Guaranteed seal and product sterility




**INDICATIONS** General soft tissue closing  
and/or ligation  
especially in Plastic surgery  
Ophthalmic surgery








**CONTRAINDICATIONS** Cardiovascular surgery  
Neurosurgery  
Microsurgery




 Round bodied taper point double armed	Metric	USP	Length	Order number
<b>1/4 circle</b>  2 x VR 8 mm 220 µm 	0.4	8/0	45 cm	<i>black</i> ON81




 Side cutting spatula	Metric	USP	Length	Order number
<b>1/2 circle</b>  HSP 5.5 mm 150 µm 	0.2	10/0	30 cm	<i>black</i> ON1081

 Side cutting spatula double armed	Metric	USP	Length	Order number
2 x HSP 5.5 mm 150 µm 	0.2	10/0	30 cm	<i>black</i> ON108
<b>1/2 circle</b>  2 x HSP 6.2 mm 150 µm 	0.2	10/0	30 cm	<i>black</i> ON102
2 x HSP 6.4 mm 220 µm 	0.2	10/0	30 cm	<i>black</i> ON106
2 x HSP 8 mm 300 µm 	0.7	6/0	45 cm	<i>black</i> ON62

 Side cutting spatula	Metric	USP	Length	Order number
<b>3/8 circle</b>  DSP 6.2 mm 150 µm 	0.2	10/0	15 cm	<i>black</i> ON1011

 Side cutting spatula double armed	Metric	USP	Length	Order number			
	2 x DSP 6 mm 220 µm		0.2	10/0	30 cm	<i>black</i>	ON103
			0.1	11/0	20 cm	<i>black</i>	ON111
			0.2	10/0	30 cm	<i>black</i>	ON101
	2 x DSP 6.2 mm 150 µm		0.3	9/0	30 cm	<i>black</i>	ON95
			0.3	9/0	45 cm	<i>black</i>	ON95-45
			0.4	8/0	30 cm	<i>black</i>	ON85
			0.4	8/0	45 cm	<i>black</i>	ON85-45
	2 x DSP 6.2 mm 200 µm		0.5	7/0	30 cm	<i>black</i>	ON75
	2 x DSP 6.4 mm 220 µm		0.2	10/0	30 cm	<i>black</i>	ON105
			0.2	10/0	38 cm	<i>black</i>	ON105-38
		0.2	10/0	45 cm	<i>black</i>	ON105-45	
2 x DSP 8 mm 300 µm		0.7	6/0	45 cm	<i>black</i>	ON61	

 Side cutting spatula double armed	Metric	USP	Length	Order number			
	2 x VSP 8 mm 350 µm		1	5/0	45 cm	<i>black</i>	ON51

 Reverse cutting	Metric	USP	Length	Order number			
	DS 12 mm		0.5	7/0	75 cm	<i>black</i>	9051512
			0.7	6/0	75 cm	<i>black</i>	9071512