


DACLON NYLON

DACLON
NYLON



Type	Monofilament
Composition	Extrusion of polyamide 6.0 or 6.6
Coating	None
Colour	Black
Tissue reaction	Minimal
Absorption	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
Presentation	Box with 12 sutures

*For quality
&
safety*

smi

www.sutures.be

DACLON

NYLON

CHARACTERISTICS


- ✓ High tensile strength, non absorbable
- ✓ Synthetic
- ✓ Monofilament with microscopic uniform diameter
- ✓ Pliability and softness
- ✓ Hermetically sealed packing

BENEFITS

- Excellent and permanent tissue support
- Minimal tissue reaction
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
- Good handling properties
- Easy and secure knotting
- Guaranteed seal and product sterility


INDICATIONS General soft tissue closing
and/or ligation
especially in Plastic surgery
Ophthalmic surgery

Round bodied taper point double armed



Metric USP Length Order number

1/4 Circle



2 x VR 8 mm
220 µm



0.4 8/0 45 cm *black* ON81

Side cutting spatula double armed



Metric USP Length Order number

1/2 Circle



2 x HSP 6,2 mm
150 µm



0.2 10/0 30 cm *black* ON102

2 x HSP 6,4 mm
220 µm



0.2 10/0 30 cm *black* ON106

2 x HSP 8 mm
300 µm



0.7 6/0 45 cm *black* ON62

Side cutting spatula double armed



Metric USP Length Order number

3/8 Circle



2 x DSP 6 mm
220 µm



0.2 10/0 30 cm *black* ON103

0.1 11/0 20 cm *black* ON111

0.2 10/0 30 cm *black* ON101

2 x DSP 6,2 mm
150 µm



0.3 9/0 30 cm *black* ON95

0.3 9/0 45 cm *black* ON95-45

0.4 8/0 30 cm *black* ON85

0.4 8/0 45 cm *black* ON85-45

2 x DSP 6,2 mm
200 µm



0.5 7/0 30 cm *black* ON75

2 x DSP 6,4 mm
220 µm



0.2 10/0 30 cm *black* ON105

0.2 10/0 38 cm *black* ON105-38

0.2 10/0 45 cm *black* ON105-45

2 x DSP 8 mm
300 µm



0.7 6/0 45 cm *black* ON61

Side cutting spatula



Metric USP Length Order number

3/8 Circle



DSP 6,2 mm
150 µm




0.2 10/0 15 cm *black* ON1011

**Side cutting spatula
 double armed**



Metric USP Length Order number

1/4 Circle



2 x VSP 8 mm
 350 µm



1 5/0 45 cm *black* ON51

Reverse cutting



Metric USP Length Order number

3/8 circle



DS 12 mm



0.5 7/0 75 cm *black* 9051512

 0.7 6/0 75 cm *black* 9071512