


DACLON

NYLON

DACLON
NYLON



Type	Monofilament
Composition	Extrusion of polyamide 6.0 or 6.6
Coating	None
Colour	Blue Black
Tissue reaction	Minimal
Absorption	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
Presentation	Box with 12 sutures

*For quality
&
safety*

smi

www.sutures.be

DACLON

NYLON

CHARACTERISTICS

- ☑ High tensile strength, non absorbable
- ☑ Synthetic
- ☑ Monofilament with microscopic uniform diameter
- ☑ Pliability and softness
- ☑ Hermetically sealed packing

BENEFITS


- Excellent and permanent tissue support
- Minimal tissue reaction
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
- Good handling properties
- Easy and secure knotting
- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation

Daclon Blue
 especially in Skin closure
 General surgery
 Plastic surgery
 Gastrointestinal surgery
 Gynaecology
 Obstetrics
 Orthopaedics

Daclon Black
 especially in Plastic surgery
 Ophthalmic surgery

Round bodied taper point double armed



Metric USP Length Order number

1/4 Circle



2 x VR 8 mm
220 µm



0.4

8/0

45 cm

black

ON81

Side cutting spatula double armed



Metric USP Length Order number

1/2 Circle



2 x HSP 6,2 mm
150 µm



0.2

10/0

30 cm

black

ON102

2 x HSP 6,4 mm
220 µm



0.2

10/0

30 cm

black

ON106

2 x HSP 8 mm
300 µm



0.7

6/0

45 cm

black

ON62

Side cutting spatula double armed



Metric USP Length Order number

3/8 Circle



2 x DSP 6 mm
220 µm



0.2

10/0

30 cm

black

ON103

0.1

11/0

20 cm

black

ON111

0.2

10/0

30 cm

black

ON101

0.3

9/0

30 cm

black

ON95

2 x DSP 6,2 mm
150 µm



0.3

9/0

45 cm

black

ON95-45

0.4

8/0

30 cm

black

ON85

0.4

8/0

45 cm

black

ON85-45

0.5

7/0

30 cm

black

ON75

2 x DSP 6,4 mm
220 µm



0.2

10/0

30 cm

black

ON105

0.2

10/0

45 cm

black

ON105-45

2 x DSP 8 mm
300 µm



0.7

6/0

45 cm

black

ON61

Side cutting spatula double armed



Metric USP Length Order number

1/4 Circle



2 x VSP 8 mm
350 µm



1

5/0

45 cm

black

ON51