


# POLYPROPYLENE

BLUE



<b>Type</b>	Monofilament
<b>Composition</b>	Polypropylene, a polymer of propylene
<b>Coating</b>	None
<b>Colour</b>	Blue
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	Non absorbable
<b>Presentation</b>	Box with 12 sutures

*For quality  
&  
safety*



**smi**

[www.sutures.be](http://www.sutures.be)

# POLYPROPYLENE

BLUE

## **CHARACTERISTICS**





- ✓ High tensile strength, non absorbable
- ✓ Resistance to repeated bending
- ✓ Synthetic
- ✓ Monofilament with microscopic uniform diameter
- ✓ Pliability and softness
- ✓ Hermetically sealed packing

**INDICATIONS** General soft tissue closing  
and/or ligation  
especially in Cardiovascular surgery  
Neurosurgery  
Ophthalmic surgery  
Microsurgery  
Plastic surgery  
Skin closure  
Orthopaedics  
Hernias (also as a mesh)  
Gastrointestinal surgery  
General surgery  
Gynaecology  
Obstetrics




## **BENEFITS**

- Excellent and permanent tissue support
- Use in cardiovascular surgery (heart prosthesis or vascular anastomosis)
- Minimal tissue reaction
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
- Good handling properties
- Easy and secure knotting
- Guaranteed seal and product sterility




 Side cutting spatula double armed	Metric	USP	Length	Order number
<b>1/2 circle</b>  2 x HSP 6.2 mm 200 µm 	0.2	10/0	30 cm	<i>blue</i> OPP106
<b>3/8 circle</b>  2 x DSP 6.4 mm 200 µm 	0.2	10/0	30 cm	<i>blue</i> OPP105
	0.3	9/0	30 cm	<i>blue</i> OPP95
 Side cutting spatula	Metric	USP	Length	Order number
<b>3/8 circle</b>  DSP 6 mm 200 µm 	0.5	7/0	45 cm	<i>blue</i> OPP71-1
 Side cutting spatula double armed	Metric	USP	Length	Order number
<b>1/4 circle</b>  2 x VSP 8 mm 300 µm 	0.7	6/0	45 cm	<i>blue</i> OPP62
 Side cutting spatula	Metric	USP	Length	Order number
<b>1/4 circle</b>  VSP 16 mm 150 µm 	0.2	10/0	10 cm	<i>blue</i> OPP103

 Side cutting spatula loop	Metric	USP	Length	Order number			
<b>1/4 circle</b> 	VSP 15.3 mm 150 µm		0.2	10/0	20 cm	<i>blue</i>	OPP113L-2
	VSP 16 mm 150 µm		0.2	10/0	20 cm	<i>blue</i>	OPP103L-2




  

 Side cutting spatula double armed	Metric	USP	Length	Order number			
<b>Straight</b> 	2 x GSP 16.15 mm 150 µm		0.2	10/0	20 cm	<i>blue</i>	OPP101
			0.3	9/0	20 cm	<i>blue</i>	OPP91

 Side cutting spatula loop	Metric	USP	Length	Order number			
<b>Straight</b> 	GSP 16.15 mm 150 µm		0.2	10/0	15 cm	<i>blue</i>	OPP110

 Side cutting spatula double armed	Metric	USP	Length	Order number			
<b>Straight</b> 	GSP 16 mm 150 µm DSP 4 mm 100 µm		0.2	10/0	10 cm	<i>blue</i>	OPP102
<b>3/8 circle</b> 