


# SILK



	<i>Silk Braided</i>	<i>Silk Virgin</i>
<b>Type</b>	Braided multifilament	Twisted multifilament
<b>Composition</b>	Braided fibres from the cocoon of the silkworm	Twisted fibres from the cocoon of the silkworm
<b>Coating</b>	Waxes Silicone	None
<b>Colour</b>	Black	Blue
<b>Tissue reaction</b>	Moderate	
<b>Absorption</b>	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.	
<b>Presentation</b>	Box with 12 sutures	

*For quality  
&  
safety*

**smi**

[www.sutures.be](http://www.sutures.be)

# SILK

## **CHARACTERISTICS**

- ☑ High tensile strength, non absorbable
- ☑ Braided or twisted structure
  
- ☑ Coated multifilament
  
- ☑ Hermetically sealed packing

## **BENEFITS**

- Good and extended tissue support for up to 3 month
  
- Excellent handling properties  
High flexibility  
High tensile strength  
Excellent knot security
  
- Soft passage through the tissues with minimal sawing,  
tissue drag and trauma  
Good knot tie down / adjustability  
Reduced capillary action
  
- Guaranteed seal and product sterility

**INDICATIONS** General soft tissue closing  
and/or ligation  
especially in General surgery  
Skin closure  
Ophthalmic surgery  
Oral surgery  
Gastrointestinal surgery  
Gynaecology  
Obstetrics

**CONTRAINDICATIONS** Urinary tract  
Biliary tract  
  
Known allergies or sensitivities to silk

# OPHTHALMIC SURGERY

**SILK**  
VIRGIN  
BLUE

## Side cutting spatula double armed



Metric

USP

Length

Order number

### 3/8 Circle



2 x DSP 6,5 mm  
200 µm



0.4

8/0

45 cm

*blue*

OS81

# OPHTHALMIC SURGERY

**SILK**  
BRAIDED  
BLACK

## Side cutting spatula double armed



Metric

USP

Length

Order number

### 3/8 Circle



2 x DSP 6,5 mm  
200 µm



0.4

8/0

45 cm

*black*

OS81B

2 x DSP 6,5 mm  
240 µm



0.5

7/0

45 cm

*black*

OS71

0.7

6/0

45 cm

*black*

OS61

2 x DSP 8 mm  
350 µm



1

5/0

30 cm

*black*

OS52

## Side cutting spatula double armed




Metric

USP

Length

Order number

### 1/4 Circle



2 x VSP 8 mm  
250 µm



0.7

6/0

45 cm

*black*

OS62

2 x VSP 10 mm  
350 µm



1

5/0

45 cm

*black*

OS51

**smi**

www.sutures.be