



SURGICRYL[®]

910


POLYGLACTINE



| | |
|------------------------|--|
| Type | Braided and coated multifilament |
| Composition | Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide |
| Coating | Poly-glycolide-co-L-lactide and calcium stearate (<1%) |
| Colour | Violet |
| Tissue reaction | Minimal |
| Absorption | The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Between 40% and 50% of tensile strength remain after 21 days. |
| Presentation | Box with 12 sutures |



*For quality
&
safety*



www.sutures.be

SURGICRYL®

910

POLYGLACTINE

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure

- Coated multifilament

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
 - High flexibility
 - High tensile strength
 - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Plastic surgery
- Ophthalmic surgery
- Gynaecology
- Obstetrics
- Orthopaedics
- Urology

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

TENSILE STRENGTH



ABSORPTION



| Side cutting spatula double armed | | Metric | USP | Length | Order number | |
|--|--------------|--------|-----|--------|---------------|------------|
| 1/2 circle  | 2 x HSP 8 mm | 0.5 | 7/0 | 45 cm | <i>violet</i> | O91007 |
| | 200 µm | 0.7 | 6/0 | 45 cm | <i>violet</i> | O910061 |
| 3/8 circle  | 2 x DSP 6 mm | 0.4 | 8/0 | 45 cm | <i>violet</i> | O91008 |
| | 200 µm | 0.5 | 7/0 | 45 cm | <i>violet</i> | O910072 |
| | 220 µm | 1 | 5/0 | 45 cm | <i>violet</i> | O910051 |
| 1/4 circle  | 2 x VSP 8 mm | 0.7 | 6/0 | 45 cm | <i>violet</i> | O91006 |
| | 300 µm | 0.7 | 6/0 | 20 cm | <i>undyed</i> | O91006U-20 |
| | | 0.7 | 6/0 | 45 cm | <i>undyed</i> | O91006U |
| | 2 x VSP 8 mm | 1 | 5/0 | 45 cm | <i>violet</i> | O91005 |
| | 350 µm | | | | | |