


# SILK



	<i>Silk Braided</i>	<i>Silk Virgin</i>
<b>Type</b>	Braided multifilament	Twisted multifilament
<b>Composition</b>	Braided fibres from the cocoon of the silkworm	Twisted fibres from the cocoon of the silkworm
<b>Coating</b>	Waxes Silicone	None
<b>Colour</b>	Black	Blue
<b>Tissue reaction</b>	Moderate	
<b>Absorption</b>	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.	
<b>Presentation</b>	Box with 12 sutures	

*For quality  
&  
safety*

**smi**

[www.sutures.be](http://www.sutures.be)

# SILK

## **CHARACTERISTICS**

- ☑ High tensile strength, non absorbable
- ☑ Braided or twisted structure
  
- ☑ Coated multifilament
  
- ☑ Hermetically sealed packing

## **BENEFITS**

- Good and extended tissue support for up to 3 month
  
- Excellent handling properties  
High flexibility  
High tensile strength  
Excellent knot security
  
- Soft passage through the tissues with minimal sawing, tissue drag and trauma  
Good knot tie down / adjustability  
Reduced capillary action
  
- Guaranteed seal and product sterility

**INDICATIONS** General soft tissue closing and/or ligation  
especially in General surgery  
Skin closure  
Ophthalmic surgery  
Oral surgery  
Gastrointestinal surgery  
Gynaecology  
Obstetrics

**CONTRAINDICATIONS** Urinary tract  
Biliary tract  
  
Known allergies or sensitivities to silk

Without needle



Metric      USP      Length      Order number

2	3/0	150 cm	<i>black</i>	8200015
3	2/0	150 cm	<i>black</i>	8300015
3.5	0	150 cm	<i>black</i>	8350015
4	1	150 cm	<i>black</i>	8400015
5	2	150 cm	<i>black</i>	8500015
6	3&4	150 cm	<i>black</i>	8600015
8	6	150 cm	<i>black</i>	8800015
2	3/0	6 x 45 cm	<i>black</i>	8200044
3	2/0	6 x 45 cm	<i>black</i>	8300044
3.5	0	6 x 45 cm	<i>black</i>	8350044
4	1	6 x 45 cm	<i>black</i>	8400044
5	2	6 x 45 cm	<i>black</i>	8500044
1.5	4/0	17 x 45 cm	<i>black</i>	8150046
2	3/0	17 x 45 cm	<i>black</i>	8200046
3	2/0	17 x 45 cm	<i>black</i>	8300046
3.5	0	17 x 45 cm	<i>black</i>	8350046
4	1	17 x 45 cm	<i>black</i>	8400046
5	2	17 x 45 cm	<i>black</i>	8500046
1.5	4/0	13 x 60 cm	<i>black</i>	8150060
2	3/0	13 x 60 cm	<i>black</i>	8200060
3	2/0	13 x 60 cm	<i>black</i>	8300060
3.5	0	13 x 60 cm	<i>black</i>	8350060
4	1	13 x 60 cm	<i>black</i>	8400060
2	3/0	12 x 75 cm	<i>black</i>	8200075
3	2/0	12 x 75 cm	<i>black</i>	8300075
3.5	0	12 x 75 cm	<i>black</i>	8350075
4	1	12 x 75 cm	<i>black</i>	8400075
5	2	12 x 75 cm	<i>black</i>	8500075

Without needle



Ligapack













Metric      USP      Length      Order number

2	3/0	250 cm	<i>black</i>	8200025
3	2/0	250 cm	<i>black</i>	8300025
3.5	0	250 cm	<i>black</i>	8350025
4	1	250 cm	<i>black</i>	8400025
5	2	250 cm	<i>black</i>	8500025
6	3&4	250 cm	<i>black</i>	8600025

Round bodied taper point 

1/2 circle




		Metric	USP	Length	Color	Order number
HR 12 mm		0.7	6/0	45 cm	black	8070112
		1	5/0	45 cm	black	8100112
		1.5	4/0	75 cm	black	8150112
HR 16 mm		0.7	6/0	45 cm	black	8070116
		1	5/0	75 cm	black	8100116
		1.5	4/0	75 cm	black	8150116
HR 20 mm		2	3/0	75 cm	black	8200116
		1	5/0	75 cm	black	8100120
		1.5	4/0	75 cm	black	8150120
HR 26 mm		2	3/0	75 cm	black	8200120
		3	2/0	75 cm	black	8300120
		3.5	0	75 cm	black	8350120
		1	5/0	75 cm	black	8100126
		1.5	4/0	75 cm	black	8150126
HR 30 mm		2	3/0	75 cm	black	8200126
		3	2/0	75 cm	black	8300126
		3.5	0	75 cm	black	8350126
		4	1	75 cm	black	8400126
		5	2	75 cm	black	8500126
		2	3/0	75 cm	black	8200130
HR 37 mm		3	2/0	75 cm	black	8300130
		3.5	0	75 cm	black	8350130
		4	1	75 cm	black	8400130
		5	2	75 cm	black	8500130
HR 40 mm		2	3/0	75 cm	black	8200137
		3	2/0	75 cm	black	8300137
		3.5	0	75 cm	black	8350137
HR 45 mm		4	1	75 cm	black	8400137
		5	2	75 cm	black	8500137
		3.5	0	75 cm	black	8350140
HR 45 mm		4	1	75 cm	black	8400140
		5	2	75 cm	black	8500140
		3.5	0	75 cm	black	8350145
HR 45 mm		4	1	75 cm	black	8400145
		5	2	75 cm	black	8500145
		3.5	0	75 cm	black	8350145

### Round bodied taper point



		Metric	USP	Length	Order number
--	--	--------	-----	--------	--------------

### 3/8 circle



DR 12 mm		1	5/0	45 cm	<i>black</i>	8100212
DR 16 mm		1.5	4/0	75 cm	<i>black</i>	8150216
		2	3/0	75 cm	<i>black</i>	8200216
DR 25 mm		2	3/0	75 cm	<i>black</i>	8200225
		3	2/0	75 cm	<i>black</i>	8300225
DR 30 mm		3	2/0	75 cm	<i>black</i>	8300230

### Round bodied taper point



		Metric	USP	Length	Order number
--	--	--------	-----	--------	--------------

### Straight

GR 40 mm		3	2/0	75 cm	<i>black</i>	8302040
----------	--	---	-----	-------	--------------	---------

### Reverse cutting



		Metric	USP	Length	Order number
--	--	--------	-----	--------	--------------

### 3/8 circle



DS 12 mm		0.5	7/0	45 cm	<i>black</i>	8051512
		0.7	6/0	45 cm	<i>black</i>	8071512
		1	5/0	45 cm	<i>black</i>	8101512
		1.5	4/0	75 cm	<i>black</i>	8151512
		2	3/0	75 cm	<i>black</i>	8201512
DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>black</i>	8071512P
DS 14 mm		0.7	6/0	45 cm	<i>black</i>	8071514
DS 16 mm		0.7	6/0	75 cm	<i>black</i>	8071516
		1	5/0	75 cm	<i>black</i>	8101516
		1.5	4/0	75 cm	<i>black</i>	8151516
		2	3/0	75 cm	<i>black</i>	8201516
		3	2/0	75 cm	<i>black</i>	8301516
DS 16 mm Easy Pass		1	5/0	50 cm	<i>black</i>	8101516P
		1.5	4/0	50 cm	<i>black</i>	8151516P
DS 16 mm black		1	5/0	75 cm	<i>black</i>	8101516VB
		1.5	4/0	75 cm	<i>black</i>	8151516VB

Reverse cutting



Metric

USP

Length

Order number

DS 19 mm



0.7	6/0	75 cm	black	8071519
1	5/0	75 cm	black	8101519
1.5	4/0	75 cm	black	8151519
2	3/0	45 cm	black	8201519-45
2	3/0	75 cm	black	8201519
3	2/0	75 cm	black	8301519
3.5	0	75 cm	black	8351519
4	1	75 cm	black	8401519

DS 19 mm  
Easy Pass



1.5	4/0	50 cm	black	8151519P
2	3/0	50 cm	black	8201519P

DS 19 mm black



1	5/0	75 cm	black	8101519VB
1.5	4/0	75 cm	black	8151519VB
2	3/0	75 cm	black	8201519VB

DS 24 mm



1	5/0	75 cm	black	8101524
1.5	4/0	75 cm	black	8151524
2	3/0	75 cm	black	8201524
3	2/0	75 cm	black	8301524
3.5	0	75 cm	black	8351524
4	1	75 cm	black	8401524

3/8 circle



DS 24 mm  
Easy Pass



2	3/0	50 cm	black	8201524P
---	-----	-------	-------	----------

DS 26 mm



2	3/0	75 cm	black	8201526
3	2/0	75 cm	black	8301526

DS 30 mm



1.5	4/0	75 cm	black	8151530
2	3/0	75 cm	black	8201530
3	2/0	75 cm	black	8301530
3.5	0	75 cm	black	8351530
4	1	75 cm	black	8401530
5	2	75 cm	black	8501530

DS 37 mm





2	3/0	75 cm	black	8201537
3	2/0	75 cm	black	8301537
3.5	0	75 cm	black	8351537
4	1	75 cm	black	8401537
5	2	75 cm	black	8501537
6	3&4	75 cm	black	8601537







DS 45 mm






3	2/0	75 cm	black	8301545
3.5	0	75 cm	black	8351545
4	1	75 cm	black	8401545
5	2	75 cm	black	8501545





Reverse cutting		Metric	USP	Length	Order number
3/8 circle	DS 60 mm 	4	1	75 cm	<i>black</i> 8401560
	DS 90 mm 	3.5	0	100 cm	<i>black</i> 8351590






Reverse cutting		Metric	USP	Length	Order number
1/2 circle	HS 12 mm 	0.7	6/0	75 cm	<i>black</i> 8071412
	HS 17 mm 	1	5/0	75 cm	<i>black</i> 8101417
		2	3/0	75 cm	<i>black</i> 8201417
	HS 20 mm 	1.5	4/0	75 cm	<i>black</i> 8151420
		2	3/0	75 cm	<i>black</i> 8201420
	HS 22 mm 	1	5/0	75 cm	<i>black</i> 8101422
		1.5	4/0	75 cm	<i>black</i> 8151422
		2	3/0	75 cm	<i>black</i> 8201422
	HS 26 mm 	3	2/0	75 cm	<i>black</i> 8301422
		1	5/0	75 cm	<i>black</i> 8101426
		1.5	4/0	75 cm	<i>black</i> 8151426
	HS 36 mm 	2	3/0	75 cm	<i>black</i> 8201426
		3	2/0	75 cm	<i>black</i> 8301426
		3.5	0	75 cm	<i>black</i> 8351426
			3.5	0	75 cm
		4	1	75 cm	<i>black</i> 8401436



Cutting		Metric	USP	Length	Order number
Straight	GS 40 mm 	3	2/0	75 cm	<i>black</i> 8302140
	GS 60 mm 	2	3/0	75 cm	<i>black</i> 8202160
		3	2/0	75 cm	<i>black</i> 8302160
		3.5	0	75 cm	<i>black</i> 8352160
		4	1	75 cm	<i>black</i> 8402160

Round bodied taper cutting point 		Metric	USP	Length	Order number	
<b>1/2 circle</b> 	HRT 18 mm 	1	5/0	75 cm	<i>black</i>	8100518
		1.5	4/0	75 cm	<i>black</i>	8150518
		2	3/0	75 cm	<i>black</i>	8200518
	HRT 27 mm 	3	2/0	75 cm	<i>black</i>	8300518
		2	3/0	75 cm	<i>black</i>	8200527
		3	2/0	75 cm	<i>black</i>	8300527

Round bodied taper cutting point 		Metric	USP	Length	Order number	
<b>3/8 circle</b> 	DRT 18 mm 	1.5	4/0	75 cm	<i>black</i>	8150418
		2	3/0	75 cm	<i>black</i>	8200418