


SURGICRYL[®]

MONOFAST

POLYGLYCAPRONE 25



Type	Monofilament
Composition	Polyglycaprone 25, a copolymer made of 75% glycolide and 25% ε-caprolactone
Coating	None
Colour	Violet / Beige - undyed
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 120 days. Approximately 60% of tensile strength remain after 7 days.
Presentation	Box with 12 sutures



*For quality
&
safety*

The logo for smi, consisting of the lowercase letters 'smi' in a bold, blue, sans-serif font with a white outline.

www.sutures.be

SURGICRYL[®]

MONOFAST

POLYGLECAPRONE 25

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Monofilament with microscopic uniform diameter

- Pliability and softness
- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down

- Good handling properties
- Easy and secure knotting

- Guaranteed seal and product sterility


INDICATIONS Superficial soft tissue and/or ligation
 especially in General surgery
 Skin closure
 Plastic surgery
 Gastrointestinal surgery
 Gynaecology
 Urology

CONTRAINDICATIONS Cardiovascular surgery
 Neurosurgery
 Ophthalmic surgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

TENSILE STRENGTH



ABSORPTION



Without needle 

Metric **USP** **Length** **Order number**

3	2/0	150 cm	<i>violet</i>	16300015
3.5	0	150 cm	<i>violet</i>	16350015
4	1	150 cm	<i>violet</i>	16400015
5	2	120 cm	<i>violet</i>	16500012
5	2	150 cm	<i>violet</i>	16500015

Without needle 

Metric **USP** **Length** **Order number**









Ligapack

















5	2	250 cm	<i>violet</i>	16500025
---	---	--------	---------------	----------








Round bodied taper point 





Metric **USP** **Length** **Order number**




1/2 Circle 	HR 13 mm 	0.7	6/0	45 cm	<i>violet</i>	16070113
		1	5/0	45 cm	<i>violet</i>	16100113
		1.5	4/0	75 cm	<i>violet</i>	16150113
	HR 17 mm 	1	5/0	45 cm	<i>violet</i>	16100117
		1.5	4/0	75 cm	<i>violet</i>	16150117
		2	3/0	75 cm	<i>violet</i>	16200117
	HR 22 mm 	1.5	4/0	75 cm	<i>violet</i>	16150122
		2	3/0	75 cm	<i>violet</i>	16200122
		3	2/0	75 cm	<i>violet</i>	16300122
	HR 26 mm 	1.5	4/0	75 cm	<i>violet</i>	16150126
		2	3/0	75 cm	<i>violet</i>	16200126
		3	2/0	75 cm	<i>violet</i>	16300126
HR 30 mm 	3.5	0	75 cm	<i>violet</i>	16350126	
	2	3/0	75 cm	<i>violet</i>	16200130	
	3	2/0	75 cm	<i>violet</i>	16300130	
HR 36 mm 	3.5	0	75 cm	<i>violet</i>	16350130	
	4	1	75 cm	<i>violet</i>	16400130	
	3	2/0	75 cm	<i>violet</i>	16300136	
HR 40 mm 	3.5	0	90 cm	<i>violet</i>	16350136	
	4	1	90 cm	<i>violet</i>	16400136	
	3	2/0	75 cm	<i>violet</i>	16300140	
		3.5	0	90 cm	<i>violet</i>	16350140
		4	1	90 cm	<i>violet</i>	16400140





Round bodied taper point 		Metric	USP	Length	Order number
1/2 Circle 	HR 48 mm 	3.5	0	90 cm	<i>violet</i> 16350148
		4	1	90 cm	<i>violet</i> 16400148
	HR 65 mm 	3	2/0	90 cm	<i>violet</i> 16300165

Reverse cutting 		Metric	USP	Length	Order number
3/8 Circle 	DS 12 mm 	0.7	6/0	45 cm	<i>undyed</i> 16071512
		0.7	6/0	75 cm	<i>undyed</i> 16071512-75
		0.7	6/0	75 cm	<i>violet</i> 16071512V
		1	5/0	45 cm	<i>undyed</i> 16101512
		1	5/0	45 cm	<i>violet</i> 16101512V
	DS 12 mm Easy Pass 	0.7	6/0	50 cm	<i>undyed</i> 16071512P
		1	5/0	50 cm	<i>undyed</i> 16101512P
		1.5	4/0	50 cm	<i>undyed</i> 16151512P
	DS 16 mm 	1	5/0	75 cm	<i>undyed</i> 16101516
		1	5/0	75 cm	<i>violet</i> 16101516V
		1.5	4/0	75 cm	<i>undyed</i> 16151516
		1.5	4/0	75 cm	<i>violet</i> 16151516V
		2	3/0	75 cm	<i>undyed</i> 16201516
	DS 16 mm Easy Pass 	1	5/0	50 cm	<i>undyed</i> 16101516P
		1.5	4/0	50 cm	<i>undyed</i> 16151516P
	DS 16 mm black 	1	5/0	75 cm	<i>violet</i> 16101516VB
	DS 19 mm 	1	5/0	45 cm	<i>undyed</i> 16101519
		1.5	4/0	75 cm	<i>undyed</i> 16151519
		1.5	4/0	75 cm	<i>violet</i> 16151519V
		2	3/0	75 cm	<i>undyed</i> 16201519
3		2/0	75 cm	<i>undyed</i> 16301519	
DS 19 mm Easy Pass 	1.5	4/0	50 cm	<i>undyed</i> 16151519P	
	2	3/0	50 cm	<i>undyed</i> 16201519P	
DS 19 mm black 	1	5/0	45 cm	<i>violet</i> 16101519VB	
	1.5	4/0	75 cm	<i>violet</i> 16151519VB	

Reverse cutting 		Metric	USP	Length		Order number
3/8 Circle 	DS 24 mm 	2	3/0	75 cm	<i>undyed</i>	16201524
		3	2/0	75 cm	<i>undyed</i>	16301524
	DS 24 mm Easy Pass 	2	3/0	75 cm	<i>undyed</i>	16201524P-75
		1.5	4/0	75 cm	<i>undyed</i>	16151526
	DS 26 mm 	2	3/0	75 cm	<i>undyed</i>	16201526
		2	3/0	75 cm	<i>violet</i>	16201526V
		3	2/0	75 cm	<i>undyed</i>	16301526
	DS 30 mm 	3	2/0	75 cm	<i>violet</i>	16301526V
		3.5	0	75 cm	<i>violet</i>	16351530V
	DS 40 mm 	3.5	0	75 cm	<i>undyed</i>	16351540
		4	1	90 cm	<i>undyed</i>	16401540

Reverse cutting 		Metric	USP	Length		Order number
1/2 Circle 	HS 36 mm 	3	2/0	75 cm	<i>undyed</i>	16301436
		3.5	0	75 cm	<i>undyed</i>	16351436
	HS 40 mm 	4	1	90 cm	<i>undyed</i>	16401440

Cutting 		Metric	USP	Length		Order number
Straight 	GS 60 mm 	2	3/0	75 cm	<i>undyed</i>	16202160
		3	2/0	75 cm	<i>undyed</i>	16302160

Blunt 		Metric	USP	Length		Order number
1/2 Circle 	HRB 40 mm 	3.5	0	90 cm	<i>violet</i>	16351840
	HRB 48 mm 	4	1	90 cm	<i>violet</i>	16401848