


CATGUT

PLAIN



Type	Twisted multifilament with a monofilament appearance
Composition	Strands of purified collagen taken from the serosal layer of selected bovines. Origin from BSE-free classified countries.
Coating	None
Colour	Ivory
Tissue reaction	Moderate
Absorption	By phagocytosis. Complete mass absorption in approximately 63 days. Approximately 50% of tensile strength remain after 7 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorbed.
Presentation	Box with 12 sutures; Cassette



*For quality
&
safety*

smi

www.sutures.be

CATGUT

PLAIN

CHARACTERISTICS

- ✓ High tensile strength
- ✓ Absorption by phagocytosis
- ✓ Twisted structure

- ✓ Uniform diameter by centerless polishing (monofilament appearance)

BENEFITS


- Excellent tissue support
- Uniform absorption
- Excellent handling properties
- High flexibility
- High tensile strength
- Excellent knot security

- Soft passage through the tissues without sawing, tissue drag and trauma
- Excellent knot tie down


INDICATIONS General soft tissue closing and/or ligation especially in
General surgery
Gastrointestinal surgery
Gynaecology
Obstetrics
Urology
Ophthalmic surgery

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery
Being absorbable not to use when an extended wound support is required.
Caution in patients with delayed wound healing (e.g. compromised immune system).
In the region of infected wounds
Known sensitivities or allergies to collagen

TENSILE STRENGTH



ABSORPTION



CATGUT

PLAIN



Without needle

Metric USP Length Order number

2	4/0	150 cm	<i>ivory</i> 1200015
3	3/0	150 cm	<i>ivory</i> 1300015
3.5	2/0	150 cm	<i>ivory</i> 1350015
4	0	150 cm	<i>ivory</i> 1400015
5	1	150 cm	<i>ivory</i> 1500015
6	2	150 cm	<i>ivory</i> 1600015
7	3	150 cm	<i>ivory</i> 1700015



Round bodied taper point

Metric USP Length Order number













1/2 circle



HR 16 mm		1.5	5/0	75 cm	<i>ivory</i> 1150116
		2	4/0	75 cm	<i>ivory</i> 1200116
HR 20 mm		2	4/0	75 cm	<i>ivory</i> 1200120
		3	3/0	75 cm	<i>ivory</i> 1300120
		3.5	2/0	75 cm	<i>ivory</i> 1350120
HR 26 mm		2	4/0	75 cm	<i>ivory</i> 1200126
		3	3/0	75 cm	<i>ivory</i> 1300126
		3.5	2/0	75 cm	<i>ivory</i> 1350126
		4	0	75 cm	<i>ivory</i> 1400126
HR 30 mm		3	3/0	75 cm	<i>ivory</i> 1300130
		3	3/0	75 cm	<i>ivory</i> 1300137
		3.5	2/0	75 cm	<i>ivory</i> 1350137
		4	0	75 cm	<i>ivory</i> 1400137
HR 37 mm		5	1	75 cm	<i>ivory</i> 1500137
		5	1	90 cm	<i>ivory</i> 1500140-90
		5	1	90 cm	<i>ivory</i> 1500140-90
HR 40 mm		4	0	75 cm	<i>ivory</i> 1400145
		5	1	75 cm	<i>ivory</i> 1500145
		6	2	75 cm	<i>ivory</i> 1600145
HR 45 mm		4	0	75 cm	<i>ivory</i> 1400145
		5	1	75 cm	<i>ivory</i> 1500145
		6	2	75 cm	<i>ivory</i> 1600145
HR 50 mm		3.5	2/0	75 cm	<i>ivory</i> 1350150

CATGUT

PLAIN

 Reverse cutting		Metric	USP	Length		Order number	
3/8 circle 	DS 12 mm		1.5	5/0	45 cm	<i>ivory</i>	1151512
			1	6/0	45 cm	<i>ivory</i>	1101516
	DS 16 mm		1.5	5/0	75 cm	<i>ivory</i>	1151516
			2	4/0	75 cm	<i>ivory</i>	1201516
			3	3/0	75 cm	<i>ivory</i>	1301516
	DS 19 mm		2	4/0	75 cm	<i>ivory</i>	1201519
			3	3/0	75 cm	<i>ivory</i>	1301519
	DS 24 mm		3.5	2/0	75 cm	<i>ivory</i>	1351519
			2	4/0	75 cm	<i>ivory</i>	1201524
	DS 24 mm		3	3/0	75 cm	<i>ivory</i>	1301524
			3.5	2/0	75 cm	<i>ivory</i>	1351524
			4	0	75 cm	<i>ivory</i>	1401524
DS 37 mm		3.5	2/0	75 cm	<i>ivory</i>	1351537	

CATGUT

PLAIN

IN CASSETTES

Metric	USP	Length	Order number
2	4/0	100 m	1020
3	3/0	100 m	1030
3.5	2/0	100 m	1035
4	0	100 m	1040
5	1	75 m	1050
6	2	50 m	1060
7	3	50 m	1070
8	4	30 m	1080
9	5	30 m	1090
10	6	30 m	1100



smi

www.sutures.be