


# DACLON NYLON

DACLON  
NYLON



<b>Type</b>	Monofilament
<b>Composition</b>	Extrusion of polyamide 6.0 or 6.6
<b>Coating</b>	None
<b>Colour</b>	Blue Black
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
<b>Presentation</b>	Box with 12 sutures; Cassette



*For quality  
&  
safety*

**smi**

[www.sutures.be](http://www.sutures.be)

# DACLON

## NYLON

### **CHARACTERISTICS**

- ☑ High tensile strength, non absorbable
- ☑ Synthetic
- ☑ Monofilament with microscopic uniform diameter
  
- ☑ Pliability and softness

### **BENEFITS**

- Excellent and permanent tissue support
- Minimal tissue reaction
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
- Good handling properties
- Easy and secure knotting

**INDICATIONS**    General soft tissue closing and/or ligation

**Daclon Blue**  
 especially in    Skin closure  
                       General surgery  
                       Plastic surgery  
                       Gastrointestinal surgery  
                       Gynaecology  
                       Obstetrics  
                       Orthopaedics

**Daclon Black**  
 especially in    Plastic surgery  
                       Ophthalmic surgery

**CONTRAINDICATIONS**    Cardiovascular surgery  
   Neurosurgery  
   Microsurgery









## Round bodied taper point


		Metric	USP	Length		Order number
HR 16 mm		1	5/0	75 cm	<i>black</i>	9100116
		1.5	4/0	75 cm	<i>blue</i>	9150116
HR 20 mm		0.7	6/0	75 cm	<i>black</i>	9070120
		1	5/0	75 cm	<i>black</i>	9100120
		1.5	4/0	75 cm	<i>blue</i>	9150120
		2	3/0	75 cm	<i>blue</i>	9200120
		3	2/0	75 cm	<i>blue</i>	9300120
HR 26 mm		1.5	4/0	75 cm	<i>blue</i>	9150126
		2	3/0	75 cm	<i>blue</i>	9200126
		3	2/0	75 cm	<i>blue</i>	9300126
		3.5	0	75 cm	<i>blue</i>	9350126
		4	1	75 cm	<i>blue</i>	9400126
HR 30 mm		3	2/0	75 cm	<i>blue</i>	9300130
		3.5	0	75 cm	<i>blue</i>	9350130
HR 37 mm		4	1	75 cm	<i>blue</i>	9400130
		3	2/0	75 cm	<i>blue</i>	9300137
HR 40 mm		3.5	0	75 cm	<i>blue</i>	9350137
		4	1	75 cm	<i>blue</i>	9400137
HR 50 mm		4	1	100 cm	<i>blue</i>	9400140
		3	2/0	75 cm	<i>blue</i>	9300140
		3.5	0	75 cm	<i>blue</i>	9350140
		4	1	75 cm	<i>blue</i>	9400140
		5	2	75 cm	<i>blue</i>	9500140



1/2 circle






 Round bodied taper point loop	Metric	USP	Length	Order number
HR 30 mm 	3.5	0	150 cm	<i>blue</i> 9350930
	4	1	150 cm	<i>blue</i> 9400930
HR 35 mm 	3.5	0	150 cm	<i>blue</i> 9350935
	4	1	150 cm	<i>blue</i> 9400935
1/2 circle  HR 40 mm 	3.5	0	150 cm	<i>blue</i> 9350940
	4	1	150 cm	<i>blue</i> 9400940
HR 50 mm 	3.5	0	150 cm	<i>blue</i> 9350950
	4	1	150 cm	<i>blue</i> 9400950

### Round bodied taper point double armed



	Metric	USP	Length	Order number
1/2 circle  2 x HR 65 mm 	3	2/0	90 cm	<i>blue</i> 9300165-2

 Blunt	Metric	USP	Length	Order number
1/2 circle  HRB 40 mm 	4	1	75 cm	<i>blue</i> 9401840



## Reverse cutting










		Metric	USP	Length		Order number
DS 12 mm		0.5	7/0	75 cm	<i>black</i>	9051512
		0.7	6/0	75 cm	<i>black</i>	9071512
		1	5/0	75 cm	<i>black</i>	9101512
		1.5	4/0	75 cm	<i>blue</i>	9151512
DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>black</i>	9071512P
		1	5/0	50 cm	<i>black</i>	9101512P
		1.5	4/0	50 cm	<i>blue</i>	9151512P
DS 16 mm		0.7	6/0	45 cm	<i>black</i>	9071516-45
		0.7	6/0	75 cm	<i>black</i>	9071516
		1	5/0	45 cm	<i>black</i>	9101516-45
		1	5/0	75 cm	<i>black</i>	9101516
		1.5	4/0	45 cm	<i>blue</i>	9151516-45
		1.5	4/0	75 cm	<i>blue</i>	9151516
DS 16 mm Easy Pass		0.7	6/0	50 cm	<i>black</i>	9071516P
		1	5/0	50 cm	<i>black</i>	9101516P
		1.5	4/0	50 cm	<i>blue</i>	9151516P
DS 16 mm black		1	5/0	75 cm	<i>black</i>	9101516VB
DS 19 mm		1	5/0	45 cm	<i>black</i>	9101519-45
		1	5/0	75 cm	<i>black</i>	9101519
		1.5	4/0	45 cm	<i>blue</i>	9151519-45
		1.5	4/0	75 cm	<i>blue</i>	9151519
		1.5	4/0	90 cm	<i>blue</i>	9151519-90
		2	3/0	45 cm	<i>blue</i>	9201519-45
		2	3/0	75 cm	<i>blue</i>	9201519
		2	3/0	90 cm	<i>blue</i>	9201519-90
		3	2/0	75 cm	<i>blue</i>	9301519
DS 19 mm Easy Pass		1.5	4/0	50 cm	<i>blue</i>	9151519P
		2	3/0	50 cm	<i>blue</i>	9201519P
		3	2/0	50 cm	<i>blue</i>	9301519P
DS 19 mm black		1	5/0	75 cm	<i>black</i>	9101519VB
		1.5	4/0	75 cm	<i>blue</i>	9151519VB

3/8 circle

















## Reverse cutting





		Metric	USP	Length		Order number
DS 24 mm		1	5/0	75 cm	<i>black</i>	9101524
		1.5	4/0	45 cm	<i>blue</i>	9151524-45
		1.5	4/0	75 cm	<i>blue</i>	9151524
		2	3/0	45 cm	<i>blue</i>	9201524-45
		2	3/0	75 cm	<i>blue</i>	9201524
		3	2/0	45 cm	<i>blue</i>	9301524-45
		3	2/0	75 cm	<i>blue</i>	9301524
		3.5	0	75 cm	<i>blue</i>	9351524
		4	1	75 cm	<i>blue</i>	9401524
		5	2	90 cm	<i>blue</i>	9501524
DS 24 mm Easy Pass		3	2/0	50 cm	<i>blue</i>	9301524P
DS 26 mm		2	3/0	45 cm	<i>blue</i>	9201526-45
		2	3/0	75 cm	<i>blue</i>	9201526
		3	2/0	45 cm	<i>blue</i>	9301526-45
		3	2/0	75 cm	<i>blue</i>	9301526
DS 30 mm		2	3/0	75 cm	<i>blue</i>	9201530
		3	2/0	75 cm	<i>blue</i>	9301530
		3.5	0	75 cm	<i>blue</i>	9351530
		4	1	75 cm	<i>blue</i>	9401530
DS 36 mm		2	3/0	75 cm	<i>blue</i>	9201536
		3	2/0	75 cm	<i>blue</i>	9301536
		3.5	0	75 cm	<i>blue</i>	9351536
		4	1	90 cm	<i>blue</i>	9401536-90
DS 40 mm		3	2/0	75 cm	<i>blue</i>	9301540
		3.5	0	75 cm	<i>blue</i>	9351540
		4	1	75 cm	<i>blue</i>	9401540
DS 45 mm		5	2	90 cm	<i>blue</i>	9501540
		3	2/0	100 cm	<i>blue</i>	9301545
		4	1	75 cm	<i>blue</i>	9401545
DS 60 mm		3	2/0	75 cm	<i>blue</i>	9301560
		3	2/0	100 cm	<i>blue</i>	9301560-100
DS 90 mm		3	2/0	100 cm	<i>blue</i>	9301590
		3.5	0	100 cm	<i>blue</i>	9351590
		4	1	100 cm	<i>blue</i>	9401590

3/8 circle



 Reverse cutting		Metric	USP	Length	Order number
	HS 17 mm 	1.5	4/0	75 cm	<i>blue</i> 9151417
		2	3/0	75 cm	<i>blue</i> 9201417
	HS 22 mm 	2	3/0	75 cm	<i>blue</i> 9201422
		HS 26 mm 	2	3/0	75 cm
	3		2/0	75 cm	<i>blue</i> 9301426
	HS 30 mm 	2	3/0	90 cm	<i>blue</i> 9201430
		3	2/0	90 cm	<i>blue</i> 9301430
	HS 36 mm 	3.5	0	75 cm	<i>blue</i> 9351430
		3.5	0	75 cm	<i>blue</i> 9351436
		4	1	75 cm	<i>blue</i> 9401436
	HS 40 mm 	5	2	75 cm	<i>blue</i> 9501436
		3.5	0	75 cm	<i>blue</i> 9351440
	HS 48 mm 	4	1	75 cm	<i>blue</i> 9401440
		4	1	100 cm	<i>blue</i> 9401448

 Cutting		Metric	USP	Length	Order number
	GS 60 mm 	2	3/0	75 cm	<i>blue</i> 9202160
		3	2/0	75 cm	<i>blue</i> 9302160
		3	2/0	75 cm	<i>blue</i> 9302175

 Cutting		Metric	USP	Length	Order number	
	NS 26 mm 	2	3/0	75 cm	<i>blue</i> 9201726	
		NS 30 mm 	3	2/0	75 cm	<i>blue</i> 9301730
			3	2/0	75 cm	<i>blue</i> 9301740

## IN CASSETTES

Metric	USP	Length	Order number
1.5	4/0	100 m	9015
2	3/0	100 m	9020
3	2/0	100 m	9030
3.5	0	100 m	9035
4	1	100 m	9040
5	2	75 m	9050
6	3&4	50 m	9060

