


# SILK



	<i>Silk Braided</i>	<i>Silk Virgin</i>
<b>Type</b>	Braided multifilament	Twisted multifilament
<b>Composition</b>	Braided fibres from the cocoon of the silkworm	Twisted fibres from the cocoon of the silkworm
<b>Coating</b>	Waxes Silicone	None
<b>Colour</b>	Black	Blue
<b>Tissue reaction</b>	Moderate	
<b>Absorption</b>	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.	
<b>Presentation</b>	Box with 12 sutures	



*For quality  
&  
safety*



[www.sutures.be](http://www.sutures.be)

Without needle



Metric USP Length Order number

2	3/0	150 cm	black	8200015
3	2/0	150 cm	black	8300015
3.5	0	150 cm	black	8350015
4	1	150 cm	black	8400015
5	2	150 cm	black	8500015
6	3&4	150 cm	black	8600015
8	6	150 cm	black	8800015
2	3/0	6 x 45 cm	black	8200044
3	2/0	6 x 45 cm	black	8300044
3.5	0	6 x 45 cm	black	8350044
4	1	6 x 45 cm	black	8400044
5	2	6 x 45 cm	black	8500044
1.5	4/0	17 x 45 cm	black	8150046
2	3/0	17 x 45 cm	black	8200046
3	2/0	17 x 45 cm	black	8300046
3.5	0	17 x 45 cm	black	8350046
4	1	17 x 45 cm	black	8400046
5	2	17 x 45 cm	black	8500046
1.5	4/0	13 x 60 cm	black	8150060
2	3/0	13 x 60 cm	black	8200060
3	2/0	13 x 60 cm	black	8300060
3.5	0	13 x 60 cm	black	8350060
4	1	13 x 60 cm	black	8400060
2	3/0	12 x 75 cm	black	8200075
3	2/0	12 x 75 cm	black	8300075
3.5	0	12 x 75 cm	black	8350075
4	1	12 x 75 cm	black	8400075
5	2	12 x 75 cm	black	8500075

Without needle



Metric USP Length Order number

2	3/0	250 cm	black	8200025
3	2/0	250 cm	black	8300025
3.5	0	250 cm	black	8350025
4	1	250 cm	black	8400025
5	2	250 cm	black	8500025
6	3&4	250 cm	black	8600025











Ligapack



## Round bodied taper point

1/2 Circle




		Metric	USP	Length	Color	Order number
HR 12 mm		0.7	6/0	45 cm	black	8070112
		1	5/0	45 cm	black	8100112
		1.5	4/0	75 cm	black	8150112
HR 16 mm		0.7	6/0	45 cm	black	8070116
		1	5/0	75 cm	black	8100116
		1.5	4/0	75 cm	black	8150116
HR 20 mm		2	3/0	75 cm	black	8200116
		1	5/0	75 cm	black	8100120
		1.5	4/0	75 cm	black	8150120
HR 26 mm		2	3/0	75 cm	black	8200120
		3	2/0	75 cm	black	8300120
		3.5	0	75 cm	black	8350120
		1	5/0	75 cm	black	8100126
		1.5	4/0	75 cm	black	8150126
HR 30 mm		2	3/0	75 cm	black	8200126
		3	2/0	75 cm	black	8300126
		3.5	0	75 cm	black	8350126
		4	1	75 cm	black	8400126
		5	2	75 cm	black	8500126
HR 37 mm		2	3/0	75 cm	black	8200130
		3	2/0	75 cm	black	8300130
		3.5	0	75 cm	black	8350130
		4	1	75 cm	black	8400130
HR 40 mm		5	2	75 cm	black	8500130
		2	3/0	75 cm	black	8200137
		3	2/0	75 cm	black	8300137
HR 45 mm		3.5	0	75 cm	black	8350137
		4	1	75 cm	black	8400137
		5	2	75 cm	black	8500137
HR 45 mm		3.5	0	75 cm	black	8350140
		4	1	75 cm	black	8400140
		5	2	75 cm	black	8500140
HR 45 mm		3.5	0	75 cm	black	8350145
		4	1	75 cm	black	8400145
		5	2	75 cm	black	8500145

### Round bodied taper point



	Metric	USP	Length	Order number
--	--------	-----	--------	--------------

### 3/8 Circle



DR 12 mm	1	5/0	45 cm	<i>black</i> 8100212
DR 16 mm	1.5	4/0	75 cm	<i>black</i> 8150216
	2	3/0	75 cm	<i>black</i> 8200216
DR 25 mm	2	3/0	75 cm	<i>black</i> 8200225
	3	2/0	75 cm	<i>black</i> 8300225
DR 30 mm	3	2/0	75 cm	<i>black</i> 8300230

### Round bodied taper point



	Metric	USP	Length	Order number
--	--------	-----	--------	--------------

### Straight

GR 40 mm	3	2/0	75 cm	<i>black</i> 8302040
----------	---	-----	-------	----------------------

### Reverse cutting



	Metric	USP	Length	Order number
--	--------	-----	--------	--------------

### 3/8 Circle



DS 12 mm	0.5	7/0	45 cm	<i>black</i> 8051512
	0.7	6/0	45 cm	<i>black</i> 8071512
	1	5/0	45 cm	<i>black</i> 8101512
	1.5	4/0	75 cm	<i>black</i> 8151512
	2	3/0	75 cm	<i>black</i> 8201512
DS 12 mm Easy Pass	0.7	6/0	50 cm	<i>black</i> 8071512P
DS 14 mm	0.7	6/0	45 cm	<i>black</i> 8071514
DS 16 mm	0.7	6/0	75 cm	<i>black</i> 8071516
	1	5/0	75 cm	<i>black</i> 8101516
	1.5	4/0	75 cm	<i>black</i> 8151516
	2	3/0	75 cm	<i>black</i> 8201516
	3	2/0	75 cm	<i>black</i> 8301516
DS 16 mm Easy Pass	1	5/0	50 cm	<i>black</i> 8101516P
	1.5	4/0	50 cm	<i>black</i> 8151516P
DS 16 mm black	1	5/0	75 cm	<i>black</i> 8101516VB
	1.5	4/0	75 cm	<i>black</i> 8151516VB

Reverse cutting



Metric

USP

Length

Order number

DS 19 mm



0.7	6/0	75 cm	black	8071519
1	5/0	75 cm	black	8101519
1.5	4/0	75 cm	black	8151519
2	3/0	45 cm	black	8201519-45
2	3/0	75 cm	black	8201519
3	2/0	75 cm	black	8301519
3.5	0	75 cm	black	8351519
4	1	75 cm	black	8401519

DS 19 mm  
Easy Pass



1.5	4/0	50 cm	black	8151519P
2	3/0	50 cm	black	8201519P

DS 19 mm black



1	5/0	75 cm	black	8101519VB
1.5	4/0	75 cm	black	8151519VB
2	3/0	75 cm	black	8201519VB

DS 24 mm



1	5/0	75 cm	black	8101524
1.5	4/0	75 cm	black	8151524
2	3/0	75 cm	black	8201524
3	2/0	75 cm	black	8301524
3.5	0	75 cm	black	8351524
4	1	75 cm	black	8401524

3/8 circle



DS 24 mm  
Easy Pass



2	3/0	50 cm	black	8201524P
---	-----	-------	-------	----------

DS 26 mm



2	3/0	75 cm	black	8201526
3	2/0	75 cm	black	8301526

DS 30 mm



1.5	4/0	75 cm	black	8151530
2	3/0	75 cm	black	8201530
3	2/0	75 cm	black	8301530
3.5	0	75 cm	black	8351530
4	1	75 cm	black	8401530
5	2	75 cm	black	8501530

DS 37 mm



2	3/0	75 cm	black	8201537
3	2/0	75 cm	black	8301537
3.5	0	75 cm	black	8351537
4	1	75 cm	black	8401537
5	2	75 cm	black	8501537
6	3&4	75 cm	black	8601537

DS 45 mm



3	2/0	75 cm	black	8301545
3.5	0	75 cm	black	8351545
4	1	75 cm	black	8401545
5	2	75 cm	black	8501545

## Reverse cutting



Metric

USP

Length

Order number

### 3/8 Circle



DS 60 mm



4

1

75 cm

*black*

8401560

DS 90 mm



3.5

0

100 cm

*black*

8351590

## Reverse cutting



Metric

USP

Length

Order number

### 1/2 Circle



HS 12 mm



0.7

6/0

75 cm

*black*

8071412

HS 17 mm



1

5/0

75 cm

*black*

8101417

1.5

4/0

75 cm

*black*

8151417

2

3/0

75 cm

*black*

8201417

HS 20 mm



1.5

4/0

75 cm

*black*

8151420

2

3/0

75 cm

*black*

8201420

HS 22 mm



1

5/0

75 cm

*black*

8101422

1.5

4/0

75 cm

*black*

8151422

2

3/0

75 cm

*black*

8201422

3

2/0

75 cm

*black*

8301422

HS 26 mm



1

5/0

75 cm

*black*

8101426

1.5

4/0

75 cm

*black*

8151426

2

3/0

75 cm

*black*

8201426

3

2/0

75 cm

*black*

8301426

3.5

0

75 cm

*black*

8351426

HS 36 mm



3.5

0

75 cm

*black*

8351436

4

1

75 cm

*black*

8401436

## Cutting



Metric

USP

Length

Order number

### Straight

GS 40 mm



3

2/0

75 cm

*black*

8302140

GS 60 mm



2

3/0

75 cm

*black*

8202160

3

2/0

75 cm

*black*

8302160

3.5

0

75 cm

*black*

8352160

4




1

75 cm



*black*

8402160





## Round bodied taper cutting point

			Metric	USP	Length	Color	Order number
<b>1/2 Circle</b> 	HRT 18 mm		1	5/0	75 cm	black	8100518
			1.5	4/0	75 cm	black	8150518
			2	3/0	75 cm	black	8200518
	HRT 27 mm		3	2/0	75 cm	black	8300518
			2	3/0	75 cm	black	8200527
			3	2/0	75 cm	black	8300527




## Round bodied taper cutting point

			Metric	USP	Length	Color	Order number
<b>3/8 Circle</b> 	DRT 18 mm		1.5	4/0	75 cm	black	8150418
			2	3/0	75 cm	black	8200418

## Side cutting spatula double armed

			Metric	USP	Length	Color	Order number		
<b>3/8 Circle</b> 	2 x DSP 6,5 mm 200 µm		0.4	8/0	45 cm	black	OS81B		
			2 x DSP 6,5 mm 240 µm		0.5	7/0	45 cm	black	OS71
					0.7	6/0	45 cm	black	OS61
			2 x DSP 8 mm 350 µm		1	5/0	30 cm	black	OS52

## Side cutting spatula double armed

			Metric	USP	Length	Color	Order number
<b>1/4 Circle</b> 	2 x VSP 8 mm 250 µm		0.7	6/0	45 cm	black	OS62
			2 x VSP 10 mm 350 µm		1	5/0	45 cm

**IN CASSETTES**

Metric	USP	Length	Order number
1.5	4/0	100 m	8015
2	3/0	100 m	8020
3	2/0	100 m	8030
3.5	0	100 m	8035
4	1	100 m	8040
5	2	100 m	8050
6	3&4	80 m	8060
7	5	60 m	8070
8	6	50 m	8080
9	7	40 m	8090
10	8	30 m	8100





Side cutting spatula  
double armed



Metric

USP

Length

Order number

3/8 Circle



2 x DSP 6,5 mm  
200 µm



0.4

8/0

45 cm

*blue*

OS81