


SUPRAMID

BLACK








| | |
|------------------------|--|
| Type | Twisted and coated multifilament Sizes 5/0 and 6/0 monofilament |
| Composition | Polyamide 6.6 strands enclosed in a polyamide 6 sheet |
| Colour | Black |
| Tissue reaction | Minimal |
| Absorption | Non absorbable, gradually encapsulated by connective tissue |
| Presentation | Box with 12 sutures |












*For quality
&
safety*





The SMI logo, consisting of the letters 'smi' in a stylized, lowercase font with a blue outline and a white fill.




www.sutures.be










| Round bodied taper point | |  | Metric | USP | Length | Order number |
|--|----------|---|----------|---|--------|----------------------|
| 1/2 Circle  | HR 16 mm |  | 1.5 | 4/0 | 75 cm | <i>black</i> 4150116 |
| | | | 2 | 3/0 | 75 cm | <i>black</i> 4200116 |
| | HR 20 mm |  | 2 | 3/0 | 45 cm | <i>black</i> 4200120 |
| | | | HR 26 mm |  | 3 | 2/0 |




| Round bodied taper point | |  | Metric | USP | Length | Order number |
|--|----------|---|--------|-----|--------|----------------------|
| Straight  | GR 22 mm |  | 1.5 | 4/0 | 50 cm | <i>black</i> 4152022 |




| Reverse cutting | |  | Metric | USP | Length | Order number |
|--|----------|---|--------|-----|--------|----------------------|
| 1/2 Circle  | HS 17 mm |  | 1.5 | 4/0 | 75 cm | <i>black</i> 4151417 |
| | | | 2 | 3/0 | 75 cm | <i>black</i> 4201417 |
| | HS 20 mm |  | 1 | 5/0 | 75 cm | <i>black</i> 4101420 |
| | | | 1.5 | 4/0 | 75 cm | <i>black</i> 4151420 |
| | HS 22 mm |  | 2 | 3/0 | 75 cm | <i>black</i> 4201420 |
| | | | 1.5 | 4/0 | 75 cm | <i>black</i> 4151422 |
| | HS 26 mm |  | 2 | 3/0 | 75 cm | <i>black</i> 4201426 |
| | | | 3 | 2/0 | 75 cm | <i>black</i> 4301426 |

| Cutting | |  | Metric | USP | Length | Order number |
|--|----------|---|--------|-----|--------|----------------------|
| Straight  | GS 20 mm |  | 1.5 | 4/0 | 50 cm | <i>black</i> 4152120 |
| | GS 25 mm |  | 1.5 | 4/0 | 60 cm | <i>black</i> 4152125 |

| Reverse cutting | |  | Metric | USP | Length | Order number |
|--|----------|---|--------|-----|--------|-------------------------|
| 3/8 circle  | DS 12 mm |  | 0.7 | 6/0 | 75 cm | <i>black</i> 4071512 |
| | | | 1 | 5/0 | 45 cm | <i>black</i> 4101512-45 |
| | | | 1 | 5/0 | 75 cm | <i>black</i> 4101512 |
| | | | 1.5 | 4/0 | 45 cm | <i>black</i> 4151512 |

| Reverse cutting  | | Metric | USP | Length | Order number |
|--|---|--------|-------|--------------------------|-------------------------|
| 3/8 circle  | DS 16 mm  | 0.7 | 6/0 | 75 cm | <i>black</i> 4071516 |
| | | 1 | 5/0 | 45 cm | <i>black</i> 4101516-45 |
| | | 1 | 5/0 | 75 cm | <i>black</i> 4101516 |
| | | 1.5 | 4/0 | 45 cm | <i>black</i> 4151516 |
| | | 1.5 | 4/0 | 75 cm | <i>black</i> 4151516-75 |
| | | 2 | 3/0 | 45 cm | <i>black</i> 4201516 |
| | | 2 | 3/0 | 75 cm | <i>black</i> 4201516-75 |
| | DS 16 mm black  | 1 | 5/0 | 75 cm | <i>black</i> 4101516VB |
| | DS 19 mm  | 0.7 | 6/0 | 75 cm | <i>black</i> 4071519 |
| | | 1 | 5/0 | 75 cm | <i>black</i> 4101519 |
| | | 1.5 | 4/0 | 45 cm | <i>black</i> 4151519 |
| | | 1.5 | 4/0 | 75 cm | <i>black</i> 4151519-75 |
| | | 2 | 3/0 | 75 cm | <i>black</i> 4201519 |
| | DS 19 mm Easy Pass  | 3 | 2/0 | 75 cm | <i>black</i> 4301519 |
| DS 19 mm black  | 1.5 | 4/0 | 75 cm | <i>black</i> 4151519P-75 | |
| | 1 | 5/0 | 75 cm | <i>black</i> 4101519VB | |
| DS 24 mm  | 1.5 | 4/0 | 75 cm | <i>black</i> 4151519VB | |
| | 2 | 3/0 | 75 cm | <i>black</i> 4201524 | |
| | 3 | 2/0 | 75 cm | <i>black</i> 4301524 | |
| DS 36 mm  | 3.5 | 0 | 75 cm | <i>black</i> 4351524 | |
| | 3 | 2/0 | 75 cm | <i>black</i> 4301536 | |
| | 3.5 | 0 | 75 cm | <i>black</i> 4351536 | |
| | | 4 | 1 | 75 cm | <i>black</i> 4401536 |

| Round bodied taper cutting point  | | Metric | USP | Length | Order number |
|--|---|--------|-----|--------|----------------------|
| 1/2 circle  | HRT 18 mm  | 1.5 | 4/0 | 75 cm | <i>black</i> 4150518 |
| | | 2 | 3/0 | 75 cm | <i>black</i> 4200518 |

| Round bodied taper cutting point  | | Metric | USP | Length | Order number |
|--|---|--------|-----|--------|----------------------|
| 3/8 circle  | DRT 18 mm  | 1.5 | 4/0 | 75 cm | <i>black</i> 4150418 |
| | | 2 | 3/0 | 75 cm | <i>black</i> 4200418 |

IN CASSETTES

| Metric | USP | Length | Order number |
|--------|-----|--------|--------------|
| 1.5 | 4/0 | 100 m | 4015 |
| 2 | 3/0 | 100 m | 4020 |
| 3 | 2/0 | 100 m | 4030 |
| 3.5 | 0 | 100 m | 4035 |
| 4 | 1 | 100 m | 4040 |
| 5 | 2 | 75 m | 4050 |
| 6 | 3&4 | 50 m | 4060 |
| 7 | 5 | 40 m | 4070 |
| 8 | 6 | 40 m | 4080 |
| 9 | 7 | 30 m | 4090 |

