


SUPRAMID

BLACK



Type	Twisted and coated multifilament Sizes 5/0 and 6/0 monofilament
Composition	Polyamide 6.6 strands enclosed in a polyamide 6 sheet
Colour	Black
Tissue reaction	Minimal
Absorption	Non absorbable, gradually encapsulated by connective tissue
Presentation	Box with 12 sutures



*For quality
&
safety*

smi

www.sutures.be

SUPRAMID

BLACK

CHARACTERISTICS

- ✓ High tensile strength, non absorbable
- ✓ Synthetic
- ✓ Twisted structure

- ✓ Multifilament in a sheath

- ✓ Pliability and softness

- ✓ Hermetically sealed packing

INDICATIONS General soft tissue closing
 and/or ligation
especially in Skin closure

BENEFITS

Excellent and permanent tissue support



















Minimal tissue reaction











Excellent handling properties
High flexibility
High tensile strength
Excellent knot security




Soft passage through the tissues without sawing,
tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action




Good handling properties
Easy and secure knotting

Guaranteed seal and product sterility

Round bodied taper point			Metric	USP	Length	Order number	
1/2 Circle 	HR 16 mm		1.5	4/0	75 cm	<i>black</i> 4150116	
			2	3/0	75 cm	<i>black</i> 4200116	
	HR 20 mm		2	3/0	45 cm	<i>black</i> 4200120	
			HR 26 mm		3	2/0	75 cm
	Round bodied taper point				Metric	USP	Length
	Straight 			GR 22 mm	1.5	4/0	50 cm
Reverse cutting			Metric	USP	Length	Order number	
1/2 Circle 	HS 17 mm		1.5	4/0	75 cm	<i>black</i> 4151417	
			2	3/0	75 cm	<i>black</i> 4201417	
	HS 20 mm		1	5/0	75 cm	<i>black</i> 4101420	
			1.5	4/0	75 cm	<i>black</i> 4151420	
	HS 22 mm		2	3/0	75 cm	<i>black</i> 4201420	
			1.5	4/0	75 cm	<i>black</i> 4151422	
	HS 26 mm		1.5	4/0	75 cm	<i>black</i> 4151426	
			2	3/0	75 cm	<i>black</i> 4201426	
			3	2/0	75 cm	<i>black</i> 4301426	
	Cutting			Metric	USP	Length	Order number
Straight 			GS 20 mm	1.5	4/0	50 cm	<i>black</i> 4152120
			GS 25 mm	1.5	4/0	60 cm	<i>black</i> 4152125

Reverse cutting		Metric	USP	Length	Order number
DS 12 mm		0.7	6/0	75 cm	<i>black</i> 4071512
		1	5/0	75 cm	<i>black</i> 4101512
		1.5	4/0	45 cm	<i>black</i> 4151512
DS 16 mm		0.7	6/0	75 cm	<i>black</i> 4071516
		1	5/0	75 cm	<i>black</i> 4101516
		1.5	4/0	45 cm	<i>black</i> 4151516
		1.5	4/0	75 cm	<i>black</i> 4151516-75
		2	3/0	45 cm	<i>black</i> 4201516
		2	3/0	75 cm	<i>black</i> 4201516-75
		DS 16 mm black		1	5/0
3/8 Circle 		0.7	6/0	75 cm	<i>black</i> 4071519
		1	5/0	75 cm	<i>black</i> 4101519
		1.5	4/0	45 cm	<i>black</i> 4151519
		1.5	4/0	75 cm	<i>black</i> 4151519-75
		2	3/0	75 cm	<i>black</i> 4201519
		3	2/0	75 cm	<i>black</i> 4301519
DS 19 mm Easy Pass		1.5	4/0	75 cm	<i>black</i> 4151519P-75
DS 19 mm black		1	5/0	75 cm	<i>black</i> 4101519VB
		1.5	4/0	75 cm	<i>black</i> 4151519VB
DS 24 mm		2	3/0	75 cm	<i>black</i> 4201524
		3	2/0	75 cm	<i>black</i> 4301524
		3.5	0	75 cm	<i>black</i> 4351524
DS 36 mm		3	2/0	75 cm	<i>black</i> 4301536
		3.5	0	75 cm	<i>black</i> 4351536
		4	1	75 cm	<i>black</i> 4401536

Round bodied taper cutting point		Metric	USP	Length	Order number
1/2 Circle 		1.5	4/0	75 cm	<i>black</i> 4150518
		2	3/0	75 cm	<i>black</i> 4200518

Round bodied taper cutting point		Metric	USP	Length	Order number
3/8 Circle 		1.5	4/0	75 cm	<i>black</i> 4150418
		2	3/0	75 cm	<i>black</i> 4200418

IN CASSETTES

Metric	USP	Length	Order number
1.5	4/0	100 m	4015
2	3/0	100 m	4020
3	2/0	100 m	4030
3.5	0	100 m	4035
4	1	100 m	4040
5	2	75 m	4050
6	3&4	50 m	4060
7	5	40 m	4070
8	6	40 m	4080
9	7	30 m	4090

