


SURGICRYL[®]

910


POLYGLACTINE



| | |
|------------------------|--|
| Type | Braided and coated multifilament |
| Composition | Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide |
| Coating | Poly-glycolide-co-L-lactide and calcium stearate (<1%) |
| Colour | Violet |
| Tissue reaction | Minimal |
| Absorption | The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Between 40% and 50% of tensile strength remain after 21 days. |
| Presentation | Box with 12 sutures |



*For quality
&
safety*



www.sutures.be

SURGICRYL®

910

POLYGLACTINE

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure

- Coated multifilament

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
 - High flexibility
 - High tensile strength
 - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Plastic surgery
- Ophthalmic surgery
- Gynaecology
- Obstetrics
- Orthopaedics
- Urology

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an extended wound support is required.


Caution in patients with delayed wound healing (e.g. compromised immune system).


TENSILE STRENGTH



ABSORPTION



| Without needle |  | Metric | USP | Length | | Order number |
|----------------|---|--------|-----|------------|---------------|--------------|
| | | 1.5 | 4/0 | 150 cm | <i>violet</i> | 15150015 |
| | | 2 | 3/0 | 150 cm | <i>violet</i> | 15200015 |
| | | 3 | 2/0 | 150 cm | <i>violet</i> | 15300015 |
| | | 3.5 | 0 | 150 cm | <i>violet</i> | 15350015 |
| | | 4 | 1 | 150 cm | <i>violet</i> | 15400015 |
| | | 5 | 2 | 150 cm | <i>violet</i> | 15500015 |
| | | 3 | 2/0 | 3 x 45 cm | <i>violet</i> | 15300045-3 |
| | | 3.5 | 0 | 3 x 45 cm | <i>violet</i> | 15350045-3 |
| | | 1.5 | 4/0 | 6 x 45 cm | <i>violet</i> | 15150045-6 |
| | | 2 | 3/0 | 6 x 45 cm | <i>violet</i> | 15200045-6 |
| | | 3 | 2/0 | 6 x 45 cm | <i>violet</i> | 15300045-6 |
| | | 3.5 | 0 | 6 x 45 cm | <i>violet</i> | 15350045-6 |
| | | 1.5 | 4/0 | 12 x 45 cm | <i>violet</i> | 15150045 |
| | | 2 | 3/0 | 12 x 45 cm | <i>violet</i> | 15200045 |
| | | 3 | 2/0 | 12 x 45 cm | <i>violet</i> | 15300045 |
| | | 3.5 | 0 | 12 x 45 cm | <i>violet</i> | 15350045 |
| | | 4 | 1 | 12 x 45 cm | <i>violet</i> | 15400045 |
| | | 5 | 2 | 12 x 45 cm | <i>violet</i> | 15500045 |
| | | 2 | 3/0 | 2 x 70 cm | <i>violet</i> | 15200070-2 |
| | | 3 | 2/0 | 2 x 70 cm | <i>violet</i> | 15300070-2 |
| | | 3.5 | 0 | 2 x 70 cm | <i>violet</i> | 15350070-2 |
| | | 4 | 1 | 2 x 70 cm | <i>violet</i> | 15400070-2 |
| | | 5 | 2 | 2 x 70 cm | <i>violet</i> | 15500070-2 |


| Without needle |  | Metric | USP | Length | | Order number |
|-----------------|---|--------|-----|--------|---------------|--------------|
| Ligapack | | 2 | 3/0 | 250 cm | <i>violet</i> | 15200025 |
| | | 3 | 2/0 | 250 cm | <i>violet</i> | 15300025 |
| | | 3.5 | 0 | 250 cm | <i>violet</i> | 15350025 |
| | | 4 | 1 | 250 cm | <i>violet</i> | 15400025 |
| | | 5 | 2 | 250 cm | <i>violet</i> | 15500025 |









Round bodied taper point









1/2 Circle











| | | Metric | USP | Length | Color | Order number |
|----------|---|--------|---------------|-------------|---------------|--------------|
| HR 13 mm | | 0.7 | 6/0 | 75 cm | <i>violet</i> | 15070113 |
| | | 1 | 5/0 | 75 cm | <i>violet</i> | 15100113 |
| | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15150113 |
| HR 17 mm | | 1 | 5/0 | 75 cm | <i>violet</i> | 15100117 |
| | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15150117 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> | 15200117 |
| HR 22 mm | | 3 | 2/0 | 75 cm | <i>violet</i> | 15300117 |
| | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15150122 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> | 15200122 |
| HR 26 mm | | 3 | 2/0 | 75 cm | <i>violet</i> | 15300122 |
| | | 4 | 1 | 75 cm | <i>violet</i> | 15400122 |
| | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15150126 |
| HR 30 mm | | 2 | 3/0 | 75 cm | <i>violet</i> | 15200126 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15300126 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15300126 |
| | | 3.5 | 0 | 75 cm | <i>violet</i> | 15350126 |
| | | 4 | 1 | 75 cm | <i>violet</i> | 15400126 |
| HR 36 mm | | 5 | 2 | 75 cm | <i>violet</i> | 15500126 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> | 15200130 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15300130 |
| | | 3.5 | 0 | 75 cm | <i>violet</i> | 15350130 |
| | | 4 | 1 | 75 cm | <i>violet</i> | 15400130 |
| HR 40 mm | | 5 | 2 | 90 cm | <i>violet</i> | 15500130 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> | 15200136 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15300136 |
| | | 3.5 | 0 | 90 cm | <i>violet</i> | 15350136 |
| | | 4 | 1 | 90 cm | <i>violet</i> | 15400136 |
| HR 45 mm | | 5 | 2 | 90 cm | <i>violet</i> | 15500136 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> | 15200140 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15300140 |
| | | 3.5 | 0 | 75 cm | <i>violet</i> | 15350140-75 |
| | | 3.5 | 0 | 90 cm | <i>violet</i> | 15350140 |
| | | 4 | 1 | 75 cm | <i>violet</i> | 15400140-75 |
| | | 4 | 1 | 90 cm | <i>violet</i> | 15400140 |
| 5 | 2 | 75 cm | <i>violet</i> | 15500140-75 | | |
| 5 | 2 | 90 cm | <i>violet</i> | 15500140 | | |












| Round bodied taper point | |  | Metric | USP | Length | Order number | | |
|--------------------------|---|---|--|-------|---------------|--------------|---------------|----------|
| HR 48 mm |  | 3 | 2/0 | 90 cm | <i>violet</i> | 15300148 | | |
| | | 3.5 | 0 | 90 cm | <i>violet</i> | 15350148 | | |
| | | 4 | 1 | 90 cm | <i>violet</i> | 15400148 | | |
| | | 5 | 2 | 90 cm | <i>violet</i> | 15500148 | | |
| HR 60 mm |  | 3.5 | 0 | 75 cm | <i>violet</i> | 15350160 | | |
| | | 4 | 1 | 75 cm | <i>violet</i> | 15400160 | | |
| 1/2 Circle |  | HR 65 mm |  | 3 | 2/0 | 90 cm | <i>violet</i> | 15300165 |
| | | | | 3.5 | 0 | 90 cm | <i>violet</i> | 15350165 |
| | | HR 80 mm |  | 4 | 1 | 100 cm | <i>violet</i> | 15400180 |

| Round bodied taper point double armed | |  | Metric | USP | Length | Order number | | |
|---------------------------------------|---|---|---|-----|--------|--------------|---------------|------------|
| 1/2 Circle |  | 2 x HR 26 mm |  | 2 | 3/0 | 75 cm | <i>violet</i> | 15200126-2 |

| Round bodied taper point loop | |  | Metric | USP | Length | Order number | | |
|-------------------------------|---|---|---|-----|--------|--------------|---------------|----------|
| 1/2 Circle |  | HR 45 mm |  | 3.5 | 0 | 150 cm | <i>violet</i> | 15350945 |

| Round bodied taper point | |  | Metric | USP | Length | | Order number |
|--|----------|---|--------|-----|--------|---------------|--------------|
| 5/8 Circle  | FR 27 mm |  | 3 | 2/0 | 75 cm | <i>violet</i> | 15300327 |
| | | | 3.5 | 0 | 75 cm | <i>violet</i> | 15350327 |

| Round bodied taper point | |  | Metric | USP | Length | | Order number |
|--|----------|---|--------|-----|--------|---------------|--------------|
| 3/8 Circle  | DR 17 mm |  | 1.5 | 4/0 | 45 cm | <i>violet</i> | 15150217 |
| | DR 22 mm |  | 2 | 3/0 | 75 cm | <i>violet</i> | 15200222 |
| | DR 30 mm |  | 4 | 1 | 75 cm | <i>violet</i> | 15400230 |

| Reverse cutting | |  | Metric | USP | Length | | Order number |
|--|---|---|--------|-------|---------------|---------------|--------------|
| 1/2 Circle  | HS 12 mm |  | 0.7 | 6/0 | 45 cm | <i>violet</i> | 15071412 |
| | HS 17 mm |  | 1 | 5/0 | 75 cm | <i>violet</i> | 15101417 |
| | | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15151417 |
| | HS 22 mm |  | 2 | 3/0 | 75 cm | <i>violet</i> | 15201422 |
| | HS 23 mm |  | 5 | 2 | 90 cm | <i>violet</i> | 15501423 |
| | HS 26 mm |  | 2 | 3/0 | 75 cm | <i>violet</i> | 15201426 |
| | | | 3 | 2/0 | 75 cm | <i>violet</i> | 15301426 |
| | | | 3.5 | 0 | 75 cm | <i>violet</i> | 15351426 |
| | | | 4 | 1 | 75 cm | <i>violet</i> | 15401426 |
| | HS 30 mm |  | 2 | 3/0 | 75 cm | <i>violet</i> | 15201430 |
| | HS 36 mm |  | 3 | 2/0 | 75 cm | <i>violet</i> | 15301436 |
| | | | 3.5 | 0 | 90 cm | <i>violet</i> | 15351436 |
| | | | 4 | 1 | 90 cm | <i>violet</i> | 15401436 |
| | | | 5 | 2 | 90 cm | <i>violet</i> | 15501436 |
| | HS 40 mm |  | 2 | 3/0 | 90 cm | <i>violet</i> | 15201440 |
| 3 | | | 2/0 | 90 cm | <i>violet</i> | 15301440 | |
| 3.5 | | | 0 | 90 cm | <i>violet</i> | 15351440 | |
| 4 | | | 1 | 90 cm | <i>violet</i> | 15401440 | |
| HS 48 mm |  | 5 | 2 | 90 cm | <i>violet</i> | 15501440 | |
| | | 3.5 | 0 | 90 cm | <i>violet</i> | 15351448 | |
| | | 4 | 1 | 90 cm | <i>violet</i> | 15401448 | |
| | | | 5 | 2 | 90 cm | <i>violet</i> | 15501448 |

Reverse cutting



Metric

USP










Length





Order number






3/8 Circle











| | | Metric | USP | Length | | Order number |
|--------------------|--|--------|-----|--------|---------------|--------------|
| DS 12 mm | | 0.7 | 6/0 | 45 cm | <i>violet</i> | 15071512 |
| | | 1 | 5/0 | 45 cm | <i>violet</i> | 15101512 |
| | | 1.5 | 4/0 | 45 cm | <i>violet</i> | 15151512 |
| DS 16 mm | | 0.7 | 6/0 | 45 cm | <i>violet</i> | 15071516 |
| | | 1 | 5/0 | 75 cm | <i>violet</i> | 15101516 |
| | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15151516 |
| | | 2 | 3/0 | 45 cm | <i>violet</i> | 15201516-45 |
| DS 16 mm Easy Pass | | 2 | 3/0 | 50 cm | <i>undyed</i> | 15201516UP |
| DS 16 mm black | | 0.7 | 6/0 | 45 cm | <i>violet</i> | 15071516VB |
| DS 19 mm | | 1 | 5/0 | 75 cm | <i>violet</i> | 15101519 |
| | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15151519 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> | 15201519 |
| DS 19 mm Easy Pass | | 1.5 | 4/0 | 50 cm | <i>undyed</i> | 15151519UP |
| DS 19 mm black | | 1 | 5/0 | 75 cm | <i>violet</i> | 15101519VB |
| | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15151519VB |
| DS 24 mm | | 1.5 | 4/0 | 75 cm | <i>violet</i> | 15151524 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> | 15201524 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15301524 |
| | | 3.5 | 0 | 75 cm | <i>violet</i> | 15351524 |
| DS 24 mm Easy Pass | | 2 | 3/0 | 50 cm | <i>undyed</i> | 15201524UP |
| DS 26 mm | | 2 | 3/0 | 75 cm | <i>violet</i> | 15201526 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15301526 |
| DS 26 mm Easy Pass | | 2 | 3/0 | 75 cm | <i>undyed</i> | 15201526UP |
| DS 30 mm | | 2 | 3/0 | 75 cm | <i>violet</i> | 15201530 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15301530 |
| | | 3.5 | 0 | 75 cm | <i>violet</i> | 15351530 |
| DS 36 mm | | 3 | 2/0 | 75 cm | <i>violet</i> | 15301536 |
| | | 3.5 | 0 | 75 cm | <i>violet</i> | 15351536 |
| | | 4 | 1 | 75 cm | <i>violet</i> | 15401536 |
| DS 40 mm | | 3.5 | 0 | 75 cm | <i>violet</i> | 15351540 |
| | | 4 | 1 | 75 cm | <i>violet</i> | 15401540 |
| DS 45 mm | | 3 | 2/0 | 75 cm | <i>violet</i> | 15301545 |

| Round bodied taper cutting point  | | Metric | USP | Length | Order number |
|--|---|--------|--------|------------------------|------------------------|
| 1/2 Circle  | HRT 13 mm  | 1.5 | 4/0 | 75 cm | <i>violet</i> 15150513 |
| | HRT 17 mm  | 1.5 | 4/0 | 75 cm | <i>violet</i> 15150517 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> 15200517 |
| | HRT 26 mm  | 1.5 | 4/0 | 75 cm | <i>violet</i> 15150526 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> 15200526 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> 15300526 |
| | HRT 30 mm  | 3 | 2/0 | 75 cm | <i>violet</i> 15300530 |
| | | 3.5 | 0 | 75 cm | <i>violet</i> 15350530 |
| | HRT 36 mm  | 4 | 1 | 75 cm | <i>violet</i> 15400530 |
| | | 3 | 2/0 | 90 cm | <i>violet</i> 15300536 |
| 3.5 | | 0 | 90 cm | <i>violet</i> 15350536 | |
| 4 | | 1 | 90 cm | <i>violet</i> 15400536 | |
| HRT 40 mm  | 3.5 | 0 | 90 cm | <i>violet</i> 15350540 | |
| | 4 | 1 | 90 cm | <i>violet</i> 15400540 | |
| | 5 | 2 | 90 cm | <i>violet</i> 15500540 | |
| HRT 48 mm  | 8 | 6 | 150 cm | <i>violet</i> 15800548 | |

| Round bodied taper cutting point  | | Metric | USP | Length | Order number |
|--|---|--------|-----|--------|------------------------|
| 3/8 Circle  | DRT 13 mm  | 1.5 | 4/0 | 75 cm | <i>violet</i> 15150413 |
| | DRT 17 mm  | 1.5 | 4/0 | 75 cm | <i>violet</i> 15150417 |
| | | 2 | 3/0 | 75 cm | <i>violet</i> 15200417 |

| Round bodied taper cutting point  | | Metric | USP | Length | Order number |
|--|---|--------|-----|--------|-----------------------------|
| 3/8 Circle  | DRT 13 mm  | 1 | 5/0 | 75 cm | <i>violet</i> 15100530-0413 |
| 1/2 Circle  | HRT 30 mm  | | | | |

| Cutting | | Metric | USP | Length | | Order number |
|--------------------------|---------------------|--------|-------|---------------|---|---|
| Straight | GS 60 mm | | | |  | |
| | | 2 | 3/0 | 75 cm | <i>violet</i> | 15202160 |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15302160 |
| Blunt | | Metric | USP | Length | | Order number |
| 1/2 Circle | HRB 30 mm | | | |  | |
| | | 3 | 2/0 | 75 cm | <i>violet</i> | 15301830 |
| | | | | | |  |
| HRB 36 mm | 3.5 | 0 | 90 cm | <i>violet</i> | 15351836 | |
| | 4 | 1 | 90 cm | <i>violet</i> | 15401836 | |
| HRB 45 mm | 5 | 2 | 90 cm | <i>violet</i> | 15501845 | |
| Blunt | | Metric | USP | Length | | Order number |
| 3/8 Circle | DRB 63 mm | | | |  | |
| | | 4 | 1 | 100 cm | <i>violet</i> | 15401963 |
| | | | | |  | |
| DRB 85 mm | 5 | 2 | 75 cm | <i>violet</i> | 15501985 | |
| | | | | |  | |
| Micro needle taper point | | Metric | USP | Length | | Order number |
| 3/8 Circle | MDR 6,4 mm 150µm | 0.4 | 8/0 | 20 cm |  | |
| | | | | |  | |
| | | | | | <i>violet</i> | M910087 |

| Side cutting spatula double armed | | | Metric | USP | Length | | Order number |
|--|--------------|---|---|-----|--------|---------------|---------------|
| 1/2 Circle  | 2 x HSP 8 mm |  | 0.5 | 7/0 | 45 cm | <i>violet</i> | O91007 |
| | 200 µm | | 0.7 | 6/0 | 45 cm | <i>violet</i> | O910061 |
| 3/8 Circle  | 2 x DSP 6 mm |  | 0.4 | 8/0 | 45 cm | <i>violet</i> | O91008 |
| | 200 µm | | 0.5 | 7/0 | 45 cm | <i>violet</i> | O910072 |
| | 220 µm | 2 x DSP 11 mm |  | 1 | 5/0 | 45 cm | <i>violet</i> |
| Side cutting spatula double armed | |  | Metric | USP | Length | | Order number |
| 1/4 circle  | 2 x VSP 8 mm |  | 0.7 | 6/0 | 45 cm | <i>violet</i> | O91006 |
| | 300 µm | | 0.7 | 6/0 | 20 cm | <i>undyed</i> | O91006U-20 |
| | | | 0.7 | 6/0 | 45 cm | <i>undyed</i> | O91006U |
| | 2 x VSP 8 mm |  | 1 | 5/0 | 45 cm | <i>violet</i> | O91005 |

IN CASSETTES

Use within 4 months after first opening

| Metric | USP | Length | Order number |
|--------|-----|--------|--------------|
| 1.5 | 4/0 | 15 m | 17151 |
| 2 | 3/0 | 15 m | 17201 |
| 2 | 3/0 | 25 m | 17202 |
| 3 | 2/0 | 15 m | 17301 |
| 3 | 2/0 | 25 m | 17302 |
| 3 | 2/0 | 50 m | 17305 |
| 3.5 | 0 | 15 m | 17351 |
| 3.5 | 0 | 25 m | 17352 |
| 3.5 | 0 | 50 m | 17355 |
| 4 | 1 | 15 m | 17401 |
| 4 | 1 | 25 m | 17402 |
| 4 | 1 | 50 m | 17405 |
| 5 | 2 | 15 m | 17501 |
| 5 | 2 | 25 m | 17502 |
| 5 | 2 | 50 m | 17505 |
| 5 | 2 | 100 m | 175010 |
| 6 | 3&4 | 15 m | 17601 |
| 6 | 3&4 | 25 m | 17602 |
| 6 | 3&4 | 50 m | 17605 |
| 7 | 5 | 15 m | 17701 |
| 7 | 5 | 25 m | 17702 |
| 7 | 5 | 50 m | 17705 |
| 8 | 6 | 15 m | 17801 |
| 8 | 6 | 25 m | 17802 |
| 8 | 6 | 50 m | 17805 |

