


SURGICRYL[®]

MONOFAST

POLYGLECAPRONE 25



Type	Monofilament
Composition	Polyglecaprone 25, a copolymer made of 75% glycolide and 25% ε-caprolactone
Coating	None
Colour	Violet / Beige - undyed
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 120 days. Approximately 60% of tensile strength remain after 7 days.
Presentation	Box with 12 sutures



*For quality
&
safety*

www.sutures.be

SURGICRYL[®]

MONOFAST

POLYGLECAPRONE 25

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Monofilament with microscopic uniform diameter

- Pliability and softness

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down

- Good handling properties
- Easy and secure knotting

- Guaranteed seal and product sterility


INDICATIONS Superficial soft tissue and/or ligation
 especially in General surgery
 Skin closure
 Plastic surgery
 Gastrointestinal surgery
 Gynaecology
 Urology

CONTRAINDICATIONS Cardiovascular surgery
 Neurosurgery
 Ophthalmic surgery


Being absorbable not to use when an extended wound support is required.


Caution in patients with delayed wound healing (e.g. compromised immune system).

TENSILE STRENGTH



ABSORPTION



Without needle		Metric	USP	Length		Order number
		3	2/0	150 cm	<i>violet</i>	16300015
		3.5	0	150 cm	<i>violet</i>	16350015
		4	1	150 cm	<i>violet</i>	16400015
		5	2	120 cm	<i>violet</i>	16500012
		5	2	150 cm	<i>violet</i>	16500015

Without needle		Metric	USP	Length		Order number
----------------	---	--------	-----	--------	--	--------------

Ligapack










5	2	250 cm	<i>violet</i>	16500025
---	---	--------	---------------	----------

Round bodied taper point		Metric	USP	Length		Order number
--------------------------	--	--------	-----	--------	--	--------------



1/2 Circle



HR 13 mm		0.7	6/0	45 cm	<i>violet</i>	16070113
		1	5/0	45 cm	<i>violet</i>	16100113
		1.5	4/0	75 cm	<i>violet</i>	16150113
HR 17 mm		1	5/0	45 cm	<i>violet</i>	16100117
		1.5	4/0	75 cm	<i>violet</i>	16150117
		2	3/0	75 cm	<i>violet</i>	16200117
HR 22 mm		1.5	4/0	75 cm	<i>violet</i>	16150122
		2	3/0	75 cm	<i>violet</i>	16200122
		3	2/0	75 cm	<i>violet</i>	16300122
HR 26 mm		1.5	4/0	75 cm	<i>violet</i>	16150126
		2	3/0	75 cm	<i>violet</i>	16200126
		3	2/0	75 cm	<i>violet</i>	16300126
HR 30 mm		3.5	0	75 cm	<i>violet</i>	16350126
		2	3/0	75 cm	<i>violet</i>	16200130
		3	2/0	75 cm	<i>violet</i>	16300130
HR 36 mm		3.5	0	75 cm	<i>violet</i>	16350130
		4	1	75 cm	<i>violet</i>	16400130
		3	2/0	75 cm	<i>violet</i>	16300136
HR 40 mm		3.5	0	90 cm	<i>violet</i>	16350136
		4	1	90 cm	<i>violet</i>	16400136
		3	2/0	75 cm	<i>violet</i>	16300140
		3.5	0	90 cm	<i>violet</i>	16350140
		4	1	90 cm	<i>violet</i>	16400140









Round bodied taper point









1/2 circle





		Metric	USP	Length		Order number
HR 48 mm 		3.5	0	90 cm	<i>violet</i>	16350148
		4	1	90 cm	<i>violet</i>	16400148
HR 65 mm 		3	2/0	90 cm	<i>violet</i>	16300165




Reverse cutting





3/8 circle

		Metric	USP	Length		Order number
DS 12 mm 		0.7	6/0	45 cm	<i>undyed</i>	16071512
		0.7	6/0	75 cm	<i>undyed</i>	16071512-75
		0.7	6/0	75 cm	<i>violet</i>	16071512V
		1	5/0	45 cm	<i>undyed</i>	16101512
DS 12 mm Easy Pass 		1	5/0	45 cm	<i>violet</i>	16101512V
		0.7	6/0	50 cm	<i>undyed</i>	16071512P
		1	5/0	50 cm	<i>undyed</i>	16101512P
DS 16 mm 		1.5	4/0	50 cm	<i>undyed</i>	16151512P
		1	5/0	75 cm	<i>undyed</i>	16101516
		1	5/0	75 cm	<i>violet</i>	16101516V
DS 16 mm Easy Pass 		1.5	4/0	75 cm	<i>undyed</i>	16151516
		1.5	4/0	75 cm	<i>violet</i>	16151516V
		2	3/0	75 cm	<i>undyed</i>	16201516
DS 16 mm black 		1	5/0	50 cm	<i>undyed</i>	16101516P
		1.5	4/0	50 cm	<i>undyed</i>	16151516P
DS 19 mm 		1	5/0	75 cm	<i>violet</i>	16101516VB
		1	5/0	45 cm	<i>undyed</i>	16101519
		1	5/0	75 cm	<i>violet</i>	16101519V
		1.5	4/0	75 cm	<i>undyed</i>	16151519
		1.5	4/0	75 cm	<i>violet</i>	16151519V
		2	3/0	75 cm	<i>undyed</i>	16201519
		2	3/0	75 cm	<i>violet</i>	16201519V
DS 19 mm Easy Pass 		3	2/0	75 cm	<i>undyed</i>	16301519
		1.5	4/0	50 cm	<i>undyed</i>	16151519P
DS 19 mm black 		2	3/0	50 cm	<i>undyed</i>	16201519P
		1	5/0	75 cm	<i>violet</i>	16101519VB
	1.5	4/0	75 cm	<i>violet</i>	16151519VB	

Reverse cutting			Metric	USP	Length	Order number
3/8 Circle 	DS 24 mm		2	3/0	75 cm	<i>undyed</i> 16201524
			3	2/0	75 cm	<i>undyed</i> 16301524
	DS 24 mm Easy Pass		2	3/0	75 cm	<i>undyed</i> 16201524P-75
			1.5	4/0	75 cm	<i>undyed</i> 16151526
	DS 26 mm		2	3/0	75 cm	<i>undyed</i> 16201526
			2	3/0	75 cm	<i>violet</i> 16201526V
			3	2/0	75 cm	<i>undyed</i> 16301526
	DS 30 mm		3	2/0	75 cm	<i>undyed</i> 16301526V
			3	2/0	75 cm	<i>violet</i> 16301526V
	DS 30 mm		3.5	0	75 cm	<i>violet</i> 16351530V
	DS 40 mm		3.5	0	75 cm	<i>undyed</i> 16351540
			4	1	90 cm	<i>undyed</i> 16401540

Reverse cutting			Metric	USP	Length	Order number
1/2 Circle 	HS 36 mm		3	2/0	75 cm	<i>undyed</i> 16301436
			3.5	0	75 cm	<i>undyed</i> 16351436
	HS 40 mm		4	1	90 cm	<i>undyed</i> 16401440

Cutting			Metric	USP	Length	Order number
Straight 	GS 60 mm		2	3/0	75 cm	<i>undyed</i> 16202160
			3	2/0	75 cm	<i>undyed</i> 16302160

Blunt			Metric	USP	Length	Order number
1/2 Circle 	HRB 40 mm		3.5	0	90 cm	<i>violet</i> 16351840
	HRB 48 mm		4	1	90 cm	<i>violet</i> 16401848

IN CASSETTES

Use within 3 months after first opening

Metric	USP	Length	Order number
1.5	4/0	15 m	16151
2	3/0	15 m	16201
2	3/0	25 m	16202
3	2/0	15 m	16301
3	2/0	25 m	16302
3.5	0	15 m	16351
3.5	0	25 m	16352
4	1	15 m	16401
4	1	25 m	16402
4	1	50 m	16405
5	2	15 m	16501
5	2	25 m	16502
5	2	50 m	16505
5	2	100 m	165010
6	3&4	15 m	16601
6	3&4	25 m	16602
6	3&4	50 m	16605

