


# SURGICRYL<sup>®</sup>

## MONOFAST

### POLYGLECAPRONE 25



<b>Type</b>	Monofilament
<b>Composition</b>	Polyglecaprone 25, a copolymer made of 75% glycolide and 25% ε-caprolactone
<b>Coating</b>	None
<b>Colour</b>	Violet / Beige - undyed
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 120 days. Approximately 60% of tensile strength remain after 7 days.
<b>Presentation</b>	Box with 12 sutures



*For quality  
&  
safety*

[www.sutures.be](http://www.sutures.be)

# SURGICRYL<sup>®</sup>

## MONOFAST

### POLYGLECAPRONE 25

#### CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Monofilament with microscopic uniform diameter
  
- Pliability and softness
  
- Hermetically sealed packing

#### BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
  
- Good handling properties
- Easy and secure knotting
  
- Guaranteed seal and product sterility


**INDICATIONS** Superficial soft tissue and/or ligation  
 especially in General surgery  
 Skin closure  
 Plastic surgery  
 Gastrointestinal surgery  
 Gynaecology  
 Urology

**CONTRAINDICATIONS** Cardiovascular surgery  
 Neurosurgery  
 Ophthalmic surgery


Being absorbable not to use when an extended wound support is required.


Caution in patients with delayed wound healing (e.g. compromised immune system).

#### TENSILE STRENGTH



#### ABSORPTION



Without needle		Metric	USP	Length		Order number
		3	2/0	150 cm	<i>violet</i>	16300015
		3.5	0	150 cm	<i>violet</i>	16350015
		4	1	150 cm	<i>violet</i>	16400015
		5	2	120 cm	<i>violet</i>	16500012
		5	2	150 cm	<i>violet</i>	16500015

Without needle		Metric	USP	Length		Order number
----------------	---	--------	-----	--------	--	--------------

**Ligapack**














5	2	250 cm	<i>violet</i>	16500025
---	---	--------	---------------	----------











Round bodied taper point		Metric	USP	Length		Order number
--------------------------	--	--------	-----	--------	--	--------------








**1/2 Circle**










<b>HR 13 mm</b>		0.7	6/0	45 cm	<i>violet</i>	16070113	
		1	5/0	45 cm	<i>violet</i>	16100113	
		1.5	4/0	75 cm	<i>violet</i>	16150113	
	<b>HR 17 mm</b>		1	5/0	45 cm	<i>violet</i>	16100117
			1.5	4/0	75 cm	<i>violet</i>	16150117
			2	3/0	75 cm	<i>violet</i>	16200117
	<b>HR 22 mm</b>		1.5	4/0	75 cm	<i>violet</i>	16150122
			2	3/0	75 cm	<i>violet</i>	16200122
			3	2/0	75 cm	<i>violet</i>	16300122
	<b>HR 26 mm</b>		1.5	4/0	75 cm	<i>violet</i>	16150126
			2	3/0	75 cm	<i>violet</i>	16200126
			3	2/0	75 cm	<i>violet</i>	16300126
<b>HR 30 mm</b>		3.5	0	75 cm	<i>violet</i>	16350126	
		2	3/0	75 cm	<i>violet</i>	16200130	
		3	2/0	75 cm	<i>violet</i>	16300130	
<b>HR 36 mm</b>		3.5	0	75 cm	<i>violet</i>	16350130	
		4	1	75 cm	<i>violet</i>	16400130	
		3	2/0	75 cm	<i>violet</i>	16300136	
<b>HR 40 mm</b>		3.5	0	90 cm	<i>violet</i>	16350136	
		4	1	90 cm	<i>violet</i>	16400136	
		3	2/0	75 cm	<i>violet</i>	16300140	
		3.5	0	90 cm	<i>violet</i>	16350140	
		4	1	90 cm	<i>violet</i>	16400140	





Round bodied taper point 		Metric	USP	Length	Order number
<b>1/2 Circle</b> 	HR 48 mm 	3.5	0	90 cm	<i>violet</i> 16350148
		4	1	90 cm	<i>violet</i> 16400148
	HR 65 mm 	3	2/0	90 cm	<i>violet</i> 16300165

Reverse cutting 		Metric	USP	Length	Order number
<b>3/8 Circle</b> 	DS 12 mm 	0.7	6/0	45 cm	<i>undyed</i> 16071512
		0.7	6/0	75 cm	<i>undyed</i> 16071512-75
		0.7	6/0	75 cm	<i>violet</i> 16071512V
		1	5/0	45 cm	<i>undyed</i> 16101512
		1	5/0	45 cm	<i>violet</i> 16101512V
	DS 12 mm Easy Pass 	0.7	6/0	50 cm	<i>undyed</i> 16071512P
		1	5/0	50 cm	<i>undyed</i> 16101512P
		1.5	4/0	50 cm	<i>undyed</i> 16151512P
	DS 16 mm 	1	5/0	75 cm	<i>undyed</i> 16101516
		1	5/0	75 cm	<i>violet</i> 16101516V
		1.5	4/0	75 cm	<i>undyed</i> 16151516
		1.5	4/0	75 cm	<i>violet</i> 16151516V
		2	3/0	75 cm	<i>undyed</i> 16201516
		2	3/0	75 cm	<i>violet</i> 16201516V
	DS 16 mm Easy Pass 	1	5/0	50 cm	<i>undyed</i> 16101516P
		1.5	4/0	50 cm	<i>undyed</i> 16151516P
	DS 16 mm black 	1	5/0	75 cm	<i>violet</i> 16101516VB
	DS 19 mm 	1	5/0	45 cm	<i>undyed</i> 16101519
		1.5	4/0	75 cm	<i>undyed</i> 16151519
		1.5	4/0	75 cm	<i>violet</i> 16151519V
2		3/0	75 cm	<i>undyed</i> 16201519	
3		2/0	75 cm	<i>undyed</i> 16301519	
DS 19 mm Easy Pass 	1.5	4/0	50 cm	<i>undyed</i> 16151519P	
	2	3/0	50 cm	<i>undyed</i> 16201519P	
DS 19 mm black 	1	5/0	45 cm	<i>violet</i> 16101519VB	
	1.5	4/0	75 cm	<i>violet</i> 16151519VB	

Reverse cutting			Metric	USP	Length	Order number
<b>3/8 Circle</b> 	DS 24 mm		2	3/0	75 cm	<i>undyed</i> 16201524
			3	2/0	75 cm	<i>undyed</i> 16301524
	DS 24 mm Easy Pass		2	3/0	75 cm	<i>undyed</i> 16201524P-75
			1.5	4/0	75 cm	<i>undyed</i> 16151526
	DS 26 mm		2	3/0	75 cm	<i>undyed</i> 16201526
			2	3/0	75 cm	<i>violet</i> 16201526V
			3	2/0	75 cm	<i>undyed</i> 16301526
	DS 30 mm		3	2/0	75 cm	<i>violet</i> 16301526V
			3.5	0	75 cm	<i>violet</i> 16351530V
	DS 40 mm		3.5	0	75 cm	<i>undyed</i> 16351540
			4	1	90 cm	<i>undyed</i> 16401540

Reverse cutting			Metric	USP	Length	Order number
<b>1/2 Circle</b> 	HS 36 mm		3	2/0	75 cm	<i>undyed</i> 16301436
			3.5	0	75 cm	<i>undyed</i> 16351436
	HS 40 mm		4	1	90 cm	<i>undyed</i> 16401440

Cutting			Metric	USP	Length	Order number
<b>Straight</b> 	GS 60 mm		2	3/0	75 cm	<i>undyed</i> 16202160
			3	2/0	75 cm	<i>undyed</i> 16302160

Blunt			Metric	USP	Length	Order number
<b>1/2 Circle</b> 	HRB 40 mm		3.5	0	90 cm	<i>violet</i> 16351840
	HRB 48 mm		4	1	90 cm	<i>violet</i> 16401848

### IN CASSETTES

Use within 2 months after first opening

Metric	USP	Length	Order number
1.5	4/0	15 m	16151
2	3/0	15 m	16201
2	3/0	25 m	16202
3	2/0	15 m	16301
3	2/0	25 m	16302
3.5	0	15 m	16351
3.5	0	25 m	16352
4	1	15 m	16401
4	1	25 m	16402
4	1	50 m	16405
5	2	15 m	16501
5	2	25 m	16502
5	2	50 m	16505
5	2	100 m	165010
6	3&4	15 m	16601
6	3&4	25 m	16602
6	3&4	50 m	16605

