


SURGICRYL[®]

PGA

POLYGLYCOLIC ACID

SURGICRYL[®]
PGA
POLYGLYCOLIC ACID



Type	Braided and coated multifilament
Composition	Polyglycolic acid
Coating	Polycaprolactone and calcium stearate (1%)
Colour	Violet / Beige - undyed
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days. Approximately 50% of tensile strength remain after 21 days.
Presentation	Box with 12 sutures



*For quality
&
safety*

smi

www.sutures.be

SURGICRYL[®]

PGA
POLYGLYCOLIC ACID

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure

- Coated multifilament

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
 - High flexibility
 - High tensile strength
 - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Gynaecology
- Obstetrics
- Plastic surgery
- Urology
- Ophthalmic surgery
- Orthopaedics

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

TENSILE STRENGTH



ABSORPTION



Without needle



Metric	USP	Length		Order number
2	3/0	150 cm	<i>violet</i>	11200015
3	2/0	150 cm	<i>violet</i>	11300015
3.5	0	150 cm	<i>violet</i>	11350015
4	1	150 cm	<i>violet</i>	11400015
5	2	150 cm	<i>violet</i>	11500015
6	3&4	150 cm	<i>violet</i>	11600015
7	5	150 cm	<i>violet</i>	11700015
8	6	150 cm	<i>violet</i>	11800015
1.5	4/0	3 x 50 cm	<i>violet</i>	11150050
2	3/0	3 x 50 cm	<i>violet</i>	11200050
3	2/0	3 x 50 cm	<i>violet</i>	11300050
3.5	0	3 x 50 cm	<i>violet</i>	11350050
4	1	3 x 50 cm	<i>violet</i>	11400050
5	2	3 x 50 cm	<i>violet</i>	11500050
2	3/0	12 x 45 cm	<i>violet</i>	11200045
3	2/0	12 x 45 cm	<i>violet</i>	11300045
3.5	0	12 x 45 cm	<i>violet</i>	11350045
4	1	12 x 45 cm	<i>violet</i>	11400045
5	2	12 x 45 cm	<i>violet</i>	11500045
1.5	4/0	6 x 45 cm	<i>violet</i>	11150045-6
2	3/0	6 x 45 cm	<i>violet</i>	11200045-6
3	2/0	6 x 45 cm	<i>violet</i>	11300045-6
3.5	0	6 x 45 cm	<i>violet</i>	11350045-6
4	1	6 x 45 cm	<i>violet</i>	11400045-6
5	2	6 x 45 cm	<i>violet</i>	11500045-6
2	3/0	10 x 60 cm	<i>violet</i>	11200060
3	2/0	10 x 60 cm	<i>violet</i>	11300060
3.5	0	10 x 60 cm	<i>violet</i>	11350060
2	3/0	5 x 70 cm	<i>violet</i>	11200070
3	2/0	5 x 70 cm	<i>violet</i>	11300070
3.5	0	5 x 70 cm	<i>violet</i>	11350070
4	1	5 x 70 cm	<i>violet</i>	11400070
5	2	5 x 70 cm	<i>violet</i>	11500070

Without needle



Ligapack



Metric	USP	Length	Order number
1.5	4/0	250 cm	<i>violet</i> 11150025
2	3/0	250 cm	<i>violet</i> 11200025
3	2/0	250 cm	<i>violet</i> 11300025
3.5	0	250 cm	<i>violet</i> 11350025
4	1	250 cm	<i>violet</i> 11400025
5	2	250 cm	<i>violet</i> 11500025
6	3&4	250 cm	<i>violet</i> 11600025














Round bodied taper point










1/2 Circle













	Metric	USP	Length	Order number
HR 13 mm	0.7	6/0	75 cm	<i>violet</i> 11070113
	1	5/0	75 cm	<i>violet</i> 11100113
	1.5	4/0	75 cm	<i>violet</i> 11150113
HR 17 mm	0.7	6/0	75 cm	<i>violet</i> 11070117
	1	5/0	75 cm	<i>violet</i> 11100117
	1.5	4/0	75 cm	<i>violet</i> 11150117
	2	3/0	75 cm	<i>violet</i> 11200117
	3	2/0	75 cm	<i>violet</i> 11300117
HR 22 mm	1.5	4/0	75 cm	<i>violet</i> 11150122
	2	3/0	75 cm	<i>violet</i> 11200122
	2	3/0	75 cm	<i>undyed</i> 11200122U
HR 26 mm	3	2/0	75 cm	<i>violet</i> 11300122
	1.5	4/0	75 cm	<i>violet</i> 11150126
	2	3/0	75 cm	<i>violet</i> 11200126
HR 30 mm	3	2/0	75 cm	<i>violet</i> 11300126
	3	2/0	75 cm	<i>undyed</i> 11300126U
	3.5	0	75 cm	<i>violet</i> 11350126
	4	1	75 cm	<i>violet</i> 11400126
	2	3/0	75 cm	<i>violet</i> 11200130
HR 36 mm	3	2/0	75 cm	<i>violet</i> 11300130
	3	2/0	75 cm	<i>undyed</i> 11300130U
	3.5	0	75 cm	<i>violet</i> 11350130
	4	1	75 cm	<i>violet</i> 11400130
	2	3/0	75 cm	<i>violet</i> 11200136
HR 36 mm	3	2/0	75 cm	<i>violet</i> 11300136
	3.5	0	90 cm	<i>violet</i> 11350136
	4	1	90 cm	<i>violet</i> 11400136

Round bodied taper point			Metric	USP	Length	Order number
1/2 Circle 	HR 40 mm		3	2/0	75 cm	<i>violet</i> 11300140
			3.5	0	90 cm	<i>violet</i> 11350140
			4	1	75 cm	<i>violet</i> 11400140-75
			4	1	90 cm	<i>violet</i> 11400140
			5	2	75 cm	<i>violet</i> 11500140-75
	HR 48 mm		5	2	90 cm	<i>violet</i> 11500140
			3	2/0	75 cm	<i>violet</i> 11300148
			3.5	0	90 cm	<i>violet</i> 11350148
			4	1	90 cm	<i>violet</i> 11400148
			5	2	90 cm	<i>violet</i> 11500148
	6	3&4	90 cm	<i>violet</i> 11600148		
Round bodied taper point			Metric	USP	Length	Order number
3/8 Circle 	DR 17 mm		1.5	4/0	75 cm	<i>violet</i> 11150217
			2	3/0	75 cm	<i>violet</i> 11200217
Round bodied taper point			Metric	USP	Length	Order number
5/8 Circle 	FR 27 mm		3	2/0	75 cm	<i>violet</i> 11300327
Round bodied taper point			Metric	USP	Length	Order number
J-needle 	JR 36 mm		3	2/0	75 cm	<i>violet</i> 11303036
			3.5	0	75 cm	<i>violet</i> 11353036
			4	1	75 cm	<i>violet</i> 11403040

Round bodied taper point		Metric	USP	Length	Order number
5/8 Circle  J-needle 	JR 36 mm FR 27 mm 	3.5	0	75 cm	<i>violet</i> 11353036-0327

Round bodied taper point loop		Metric	USP	Length	Order number
1/2 Circle 	HR 50 mm 	4	1	150 cm	<i>violet</i> 11400950

Round bodied taper cutting point		Metric	USP	Length	Order number
1/2 Circle 	HRT 26 mm 	3	2/0	75 cm	<i>violet</i> 11300526
	HRT 30 mm 	3	2/0	75 cm	<i>violet</i> 11300530
		3.5	0	75 cm	<i>violet</i> 11350530
		4	1	75 cm	<i>violet</i> 11400530
	HRT 36 mm 	3	2/0	90 cm	<i>violet</i> 11300536
		3.5	0	90 cm	<i>violet</i> 11350536
		4	1	90 cm	<i>violet</i> 11400536
	HRT 40 mm 	3.5	0	90 cm	<i>violet</i> 11350540
		4	1	90 cm	<i>violet</i> 11400540
	HRT 48 mm 	5	2	90 cm	<i>violet</i> 11500540
HRT 48 mm	4	1	90 cm	<i>violet</i> 11400548	

Round bodied taper cutting point		Metric	USP	Length	Order number
J-needle 	JRT 31 mm 	4	1	75 cm	<i>violet</i> 11403131




Reverse cutting










3/8 Circle



		Metric	USP	Length		Order number
DS 12 mm		0.7	6/0	45 cm	<i>violet</i>	11071512
		0.7	6/0	45 cm	<i>undyed</i>	11071512U
		1	5/0	45 cm	<i>violet</i>	11101512
		1.5	4/0	75 cm	<i>violet</i>	11151512
DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>undyed</i>	11071512UP
		1	5/0	50 cm	<i>undyed</i>	11101512UP
		1.5	4/0	50 cm	<i>undyed</i>	11151512UP
DS 16 mm		0.7	6/0	75 cm	<i>violet</i>	11071516
		1	5/0	75 cm	<i>violet</i>	11101516
		1	5/0	75 cm	<i>undyed</i>	11101516U
		1.5	4/0	75 cm	<i>violet</i>	11151516
		1.5	4/0	75 cm	<i>undyed</i>	11151516U
		2	3/0	75 cm	<i>violet</i>	11201516
DS 16 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i>	11151516UP
DS 19 mm		1	5/0	75 cm	<i>violet</i>	11101519
		1.5	4/0	75 cm	<i>violet</i>	11151519
		1.5	4/0	75 cm	<i>undyed</i>	11151519U
		2	3/0	75 cm	<i>violet</i>	11201519
		2	3/0	75 cm	<i>undyed</i>	11201519U
		3	2/0	75 cm	<i>violet</i>	11301519
DS 19 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i>	11151519UP
		2	3/0	50 cm	<i>undyed</i>	11201519UP
DS 24 mm		2	3/0	75 cm	<i>violet</i>	11201524
		2	3/0	75 cm	<i>undyed</i>	11201524U
		3	2/0	75 cm	<i>violet</i>	11301524
		3	2/0	75 cm	<i>undyed</i>	11301524U
		3.5	0	75 cm	<i>violet</i>	11351524
DS 24 mm Easy Pass		2	3/0	50 cm	<i>undyed</i>	11201524UP
		3	2/0	50 cm	<i>undyed</i>	11301524UP
DS 30 mm		1.5	4/0	75 cm	<i>violet</i>	11151530
		2	3/0	75 cm	<i>violet</i>	11201530
		3	2/0	75 cm	<i>violet</i>	11301530
		3	2/0	75 cm	<i>undyed</i>	11301530U
		3.5	0	75 cm	<i>violet</i>	11351530
		4	1	75 cm	<i>violet</i>	11401530
		5	2	75 cm	<i>violet</i>	11501530

Reverse cutting			Metric	USP	Length	Order number
3/8 Circle 	DS 36 mm					
			3	2/0	75 cm	<i>violet</i> 11301536
			3.5	0	75 cm	<i>violet</i> 11351536
			4	1	75 cm	<i>violet</i> 11401536
			5	2	75 cm	<i>violet</i> 11501536
	DS 40 mm					
			3	2/0	75 cm	<i>violet</i> 11301540
			3.5	0	75 cm	<i>violet</i> 11351540
			4	1	75 cm	<i>violet</i> 11401540
			5	2	75 cm	<i>violet</i> 11501540

Reverse cutting			Metric	USP	Length	Order number
1/2 Circle 	HS 22 mm					
			2	3/0	75 cm	<i>violet</i> 11201422
			2	3/0	75 cm	<i>undyed</i> 11201422U
			3	2/0	75 cm	<i>violet</i> 11301422
	HS 26 mm					
			3	2/0	75 cm	<i>undyed</i> 11301422U
			1.5	4/0	75 cm	<i>violet</i> 11151426
			2	3/0	75 cm	<i>violet</i> 11201426
	HS 30 mm					
			3	2/0	75 cm	<i>violet</i> 11301426
			3	2/0	75 cm	<i>violet</i> 11301426
	HS 36 mm					
			3.5	0	75 cm	<i>violet</i> 11351426
			4	1	90 cm	<i>violet</i> 11401436
			5	2	90 cm	<i>violet</i> 11501436
	HS 40 mm					
			3	2/0	90 cm	<i>violet</i> 11301440
			3.5	0	90 cm	<i>violet</i> 11351440
			4	1	90 cm	<i>violet</i> 11401440
	HS 48 mm					
5			2	90 cm	<i>violet</i> 11501440	
4			1	90 cm	<i>violet</i> 11401448	
5			2	90 cm	<i>violet</i> 11501448	
			6	3&4	90 cm	<i>violet</i> 11601448

Cutting		Metric	USP	Length		Order number	
Straight	GS 60 mm						
			2	3/0	75 cm	<i>violet</i>	11202160
			2	3/0	75 cm	<i>undyed</i>	11202160U
		3	2/0	75 cm	<i>violet</i>	11302160	

Cutting		Metric	USP	Length		Order number
3/8 Circle	NS 16 mm	1	5/0	45 cm	<i>violet</i>	11101716

Side cutting spatula double armed		Metric	USP	Length		Order number
1/2 Circle	2 x HSP 8 mm 200 µm	0.5	7/0	45 cm	<i>violet</i>	OSU07

Side cutting spatula double armed		Metric	USP	Length		Order number
3/8 Circle	2 x DSP 6 mm 200 µm	0.4	8/0	45 cm	<i>violet</i>	OSU08
	2 x DSP 6,2 mm 150 µm	0.2	10/0	30 cm	<i>violet</i>	OSU10
	2 x DSP 6,4 mm 200 µm	0.3	9/0	30 cm	<i>violet</i>	OSU09
	2 x DSP 8 mm 200 µm	0.5	7/0	45 cm	<i>violet</i>	OSU071

Side cutting spatula		Metric	USP	Length		Order number
3/8 Circle	DSP 6,2 mm 150 µm	0.2	10/0	10 cm	<i>violet</i>	OSU101

Side cutting spatula double armed		Metric	USP	Length		Order number
1/4 Circle	2 x VSP 8 mm 300 µm	0.7	6/0	45 cm	<i>violet</i>	OSU06
	350 µm	1	5/0	45 cm	<i>violet</i>	OSU05

IN CASSETTES

Use within 4 months after first opening

Metric	USP	Length	Order number
1.5	4/0	15 m	11151
2	3/0	15 m	11201
2	3/0	25 m	11202
3	2/0	15 m	11301
3	2/0	25 m	11302
3	2/0	50 m	11305
3.5	0	15 m	11351
3.5	0	25 m	11352
3.5	0	50 m	11355
4	1	15 m	11401
4	1	25 m	11402
4	1	50 m	11405
5	2	15 m	11501
5	2	25 m	11502
5	2	50 m	11505
5	2	100 m	115010
6	3&4	15 m	11601
6	3&4	25 m	11602
6	3&4	50 m	11605
6	3&4	100 m	116010
7	5	15 m	11701
7	5	25 m	11702
7	5	50 m	11705
8	6	15 m	11801
8	6	25 m	11802
8	6	50 m	11805



with storage solution

Metric	USP	Length	Order number
1.5	4/0	15 m	15151
2	3/0	15 m	15201
2	3/0	25 m	15202
3	2/0	15 m	15301
3	2/0	25 m	15302
3	2/0	50 m	15305
3.5	0	15 m	15351
3.5	0	25 m	15352
3.5	0	50 m	15355
4	1	15 m	15401
4	1	25 m	15402
4	1	50 m	15405
5	2	15 m	15501

Metric	USP	Length	Order number
5	2	25 m	15502
5	2	50 m	15505
5	2	100 m	155010
6	3&4	15 m	15601
6	3&4	25 m	15602
6	3&4	50 m	15605
6	3&4	100 m	156010
7	5	15 m	15701
7	5	25 m	15702
7	5	50 m	15705
8	6	15 m	15801
8	6	25 m	15802
8	6	50 m	15805

Use within 4 months after first opening