

EPISIO SET



Product

This special set contains all sutures for a episiotomy procedures in one single foil pack.

Benefits

- Easy and quick identification of material
- All in one needle suture package

Features

- Supple passage through tissue
- Fast absorption
- High tensile strength

Technical description

Surgicryl Rapid: synthetic absorbable Polyglycolic Acid - braided and coated multifilament
12 pc in a box


Reverse cutting

Round bodied taper cutting point



Order no

| | | | Order no |
|------------|-----------|--------------------------------|---------------|
| 3/8 Circle | DS 26 mm | 1 x USP 2/0 – DS 26 mm – 75 cm | S-1430-1435 |
| | HRT 40 mm | 2 x USP 0 – HRT 40 mm – 90 cm | |
| 1/2 Circle | DS 24 mm | 1 x USP 2/0 – DS 24 mm – 75 cm | S-1430-1435-1 |
| | HRT 40 mm | 1 x USP 0 – HRT 40 mm – 90 cm | |



For quality
&
safety

smi
www.sutures.be