

# FLEX-BANDAGE



## Product

A strong bandage for use on large and small animals when maximum durability and support is required

## Use

For multiple applications in medical bandaging, ideal for first aid, running bandages, tail wrap or to wrap packed hooves

## Product highlights

- ✓ Non-slip bandage adheres to itself – stretchable - lightweight
- ✓ Maintains original size and does not shrink
- ✓ Easily tears by hand to required length
- ✓ No clips, pins or fasteners needed
- ✓ Sweat and water resistant
- ✓ Abrasion resistant to help prevent rundown

## Attributes

- ✓ Non-irritating, won't stick to skin and hair
- ✓ Does not absorb moisture or harbour bacteria growth
- ✓ Excellent porosity and allows the skin to breathe
- ✓ Excellent pressure bandage without cutting off circulation
- ✓ Contains latex

For quick removal bandage can be unwrapped or cut with scissors

## Option no bit no chew (with bitter substances)

Prevents biting and chewing of the bandages

## Technical description

Measurements: width from 2,5 cm to 15 cm length 4,5m/roll

12 colours: yellow, neon orange, neon pink, neon green, blue, black, green, white, red, smiley blue, smiley yellow, paw white

## Packing


single colour pack  
multicolour pack



**For quality  
&  
safety**

**smi**  
www.sutures.be

# FLEX-BANDAGE



2.5 cm - 4.5 m	Unit	Order no.
red	60 pc.	<b>ZBAR25</b>
green	60 pc.	<b>ZBAG25</b>
blue	60 pc.	<b>ZBAB25</b>
yellow	60 pc.	<b>ZBAY25</b>
yellow smiley	60 pc.	<b>ZBAY25-S</b>

7.5 cm - 4.5 m	Unit	Order no.
white	24 pc.	<b>ZBAW75</b>
red	24 pc.	<b>ZBAR75</b>
green	24 pc.	<b>ZBAG75</b>
blue	24 pc.	<b>ZBAB75</b>
yellow	24 pc.	<b>ZBAY75</b>
neon orange	24 pc.	<b>ZBANO75</b>
neon pink	24 pc.	<b>ZBANP75</b>
neon green	24 pc.	<b>ZBANG75</b>
blue smiley	24 pc.	<b>ZBAB75-S</b>
yellow smiley	24 pc.	<b>ZBAY75-S</b>
white paw	24 pc.	<b>ZBAW75-P</b>
different colours	24 pc.	<b>ZBAM75</b>
fun box with smiley and paw	24 pc.	<b>ZBAM75-S</b>
no bit no chew	24 pc.	<b>ZBABT75</b>

15 cm - 4.5 m	Unit	Order no.
red	12 pc.	<b>ZBAR15</b>
green	12 pc.	<b>ZBAG15</b>

5 cm - 4.5 m	Unit	Order no.
white	36 pc.	<b>ZBAW5</b>
red	36 pc.	<b>ZBAR5</b>
green	36 pc.	<b>ZBAG5</b>
blue	36 pc.	<b>ZBAB5</b>
yellow	36 pc.	<b>ZBAY5</b>
neon orange	36 pc.	<b>ZBANO5</b>
neon pink	36 pc.	<b>ZBANP5</b>
neon green	36 pc.	<b>ZBANG5</b>
blue smiley	36 pc.	<b>ZBAB5-S</b>
yellow smiley	36 pc.	<b>ZBAY5-S</b>
white paw	36 pc.	<b>ZBAW5-P</b>
different colours	36 pc.	<b>ZBAM5</b>
fun box with smiley and paw	36 pc.	<b>ZBAM5-S</b>
no bit no chew	36 pc.	<b>ZBABT5</b>

10 cm - 4.5 m	Unit	Order no.
white	18 pc.	<b>ZBAW10</b>
red	18 pc.	<b>ZBAR10</b>
green	18 pc.	<b>ZBAG10</b>
blue	18 pc.	<b>ZBAB10</b>
yellow	18 pc.	<b>ZBAY10</b>
black	18 pc.	<b>ZBABL10</b>
neon orange	18 pc.	<b>ZBANO10</b>
neon pink	18 pc.	<b>ZBANP10</b>
neon green	18 pc.	<b>ZBANG10</b>
blue smiley	18 pc.	<b>ZBAB10-S</b>
yellow smiley	18 pc.	<b>ZBAY10-S</b>
white paw	18 pc.	<b>ZBAW10-P</b>
different colours	18 pc.	<b>ZBAM10</b>
fun box with smiley and paw	18 pc.	<b>ZBAM10-S</b>
no bit no chew	18 pc.	<b>ZBABT10</b>