




MATERIAL FOR REPAIRING THE CRUCIATE LIGAMENTS

Product

With this system you can unassisted and without much effort perform a solid and long lasting fixation of the cruciate repair.




Needle



Nylon surgical
sutures



Crimp tube



Crimp forceps



Compound crimp high force

Attributes

- ☒ Stainless steel material
- ☒ Individually packed
- ☒ Useable for all dog weights
- ☒ Ethylene oxide sterilization of surgical suture

*For quality
&
safety*

smi
www.sutures.be

CRUCIATE NEEDLE

		Unit	Order no
size 2 large	length 80 mm	1 dz.	ZC0000002
size 2 large	length 80 mm	6 pc.	ZC0000002-6
size 3 medium	length 62 mm	1 dz.	ZC0000003
size 3 medium	length 62 mm	6 pc.	ZC0000003-6
size 4 small	length 45 mm	1 dz.	ZC0000004
size 4 small	length 45 mm	6 pc.	ZC0000004-6
size 5 very small	length 30 mm	1 dz.	ZC0000005
size 5 very small	length 30 mm	6 pc.	ZC0000005-6

CRUCIATE REPAIR CRIMP TUBE

	Unit	Order no
10 mm (SS 316)	1 dz.	ZC00010
12 mm (SS 316)	1 dz.	ZC00012
14 mm (SS 316)	1 dz.	ZC00014

NYLON CRUCIATE SURGICAL SUTURES

	Unit	Order no
CCL 50 lb x 500 mm dogs < 15 kg	1 pc.	ZC050
CCL 80 lb x 800 mm dogs < 20 kg	1 pc.	ZC080
CCL 100 lb x 800 mm dogs > 20 kg dogs > 40 Kg = 2 x 100 Lb	1 pc.	ZC100

NYLON CRUCIATE SURGICAL SUTURES WITH 1 CRIMP

	Unit	Order no
CCL 50 lb x 500 mm + 10 mm crimp	1 pc.	ZC05010
CCL 80 lb x 800 mm + 12 mm crimp	1 pc.	ZC08012
CCL 100 lb x 800 mm + 12 mm crimp	1 pc.	ZC10012

NYLON CRUCIATE SURGICAL SUTURES WITH 3 CRIMPS

	Unit	Order no
CCL 50 lb x 500 mm + 3 x 10 mm crimp	1 pc.	ZC05010-3
CCL 80 lb x 800 mm + 3 x 12 mm crimp	1 pc.	ZC08012-3
CCL 100 lb x 800 mm + 3 x 12 mm crimp	1 pc.	ZC10012-3

NYLON CRUCIATE SURGICAL SUTURES WITH NEEDLE + CRIMP

	Unit	Order no
CCL 50 lb x 500 mm + very small 30 mm needle + 10 mm crimp	1 pc.	ZC0501005
CCL 80 lb x 800 mm + small 45 mm needle + 12 mm crimp	1 pc.	ZC0801204
CCL 100 lb x 800 mm + medium 62 mm needle + 12 mm crimp	1 pc.	ZC1001203

NYLON CRUCIATE SURGICAL SUTURES LOOP (DOUBLE LINE)

	Unit	Order no
CCL 50 lb x 1000 mm + very small 30 mm needle loop	1 pc.	ZC0500005L
CCL 80 lb x 1600 mm + small 45 mm needle loop	1 pc.	ZC0800004L
CCL 80 lb x 1600 mm + medium 62 mm needle loop	1 pc.	ZC0800003L
CCL 100 lb x 1600 mm + medium 62 mm needle loop	1 pc.	ZC1000003L

NYLON CRUCIATE SURGICAL SUTURES WITH 2 CRIMPS LOOP (DOUBLE LINE)

	Unit	Order no
CCL 50 lb x 1000 mm + very small 30 mm needle + 2 x 10 mm crimps loop	1 pc.	ZC0501005L
CCL 80 lb x 1600 mm + small 45 mm needle + 2 x 12 mm crimps loop	1 pc.	ZC0801204L
CCL 100 lb x 1600 mm + medium 62 mm needle + 2 x 12 mm crimps loop	1 pc.	ZC1001203L

CCL INSTRUMENTS

	Unit	Order no
Crimp forceps 18 cm	1 pc.	ZC24
Compound crimp high force 23 cm	1 pc.	ZC25