

# POLYPROPYLENE

WHITE



## SURGICAL MESH

### Product

Surgical mesh is a medical device that is used to provide additional support to weakened or damaged tissue.

### Indications

Reinforcement of the abdominal wall hernia / eventration via celioscopy or laparotomy

### Attributes

- ✓ No shrinkage
- ✓ No coating
- ✓ Ethylene Oxide sterilisation
- ✓ Non-absorbable

### Benefits

- ✓ Extreme dimensional stability
- ✓ Very thin mesh structure
- ✓ Ideal porosity – well tolerated
- ✓ Excellent transparency
- ✓ Good handling - easy to cut

### Technical composition

Monofilament polypropylene knitted into an elastic, durable mesh

### Square & Rectangular Mesh

|                 | Size           | Order no        |
|-----------------|----------------|-----------------|
| standard mesh   | 5 cm x 10 cm   | <b>ZMS0510</b>  |
| 5 pieces / pack | 6 cm x 11 cm   | <b>ZMS0611</b>  |
|                 | 7,5 cm x 15 cm | <b>ZMS07515</b> |
|                 | 10 cm x 15 cm  | <b>ZMS1015</b>  |
|                 | 15 cm x 15 cm  | <b>ZMS1515</b>  |
|                 | 30 cm x 30 cm  | <b>ZMS3030</b>  |



*For quality  
&  
safety*

**smi**  
[www.sutures.be](http://www.sutures.be)