


SURGICAL MESH



Product

Surgical mesh is a medical device that is used to provide additional support to weakened or damaged tissue.

Indications

Reinforcement of the abdominal wall hernia / eventration via celioscopy or laparotomy

Attributes

- No shrinkage
- No coating
- Ethylene Oxide sterilisation
- Non-absorbable

Benefits

- Extreme dimensional stability
- Very thin mesh structure
- Ideal porosity – well tolerated
- Excellent transparency
- Good handling - easy to cut

Technical composition

Monofilament polypropylene knitted into an elastic, durable mesh

Square & Rectangular Mesh

standard mesh
5 pieces / pack

Size

Order no

5 cm x 10 cm	ZMS0510
6 cm x 11 cm	ZMS0611
7,5 cm x 15 cm	ZMS07515
10 cm x 15 cm	ZMS1015
15 cm x 15 cm	ZMS1515
30 cm x 30 cm	ZMS3030



**For quality
&
safety**



www.sutures.be