


POLYPROPYLENE

BLUE



Type	Monofilament
Composition	Polypropylene, a polymer of propylene
Coating	None
Colour	Blue
Tissue reaction	Minimal
Absorption	Non absorbable
Presentation	Box with 12 sutures

*For quality
&
safety*

smi

www.sutures.be

POLYPROPYLENE

BLUE

CHARACTERISTICS

- ✓ High tensile strength, non absorbable
- ✓ Synthetic
- ✓ Monofilament with microscopic uniform diameter

- ✓ Pliability and softness

- ✓ Hermetically sealed packing

BENEFITS

- Excellent and permanent tissue support
- Minimal tissue reaction
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down

- Good handling properties
- Easy and secure knotting

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation
especially in Ophthalmic surgery
Plastic surgery
Skin closure
Orthopaedics
Hernias (also as a mesh)
Gastrointestinal surgery
General surgery
Gynaecology
Obstetrics









CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery
Microsurgery



Round bodied taper point

1/2 circle




			Metric	USP	Length		Order number
HR 10 mm			0.5	7/0	45 cm	<i>blue</i>	5050110
			0.7	6/0	75 cm	<i>blue</i>	5070110
HR 15 mm			0.7	6/0	75 cm	<i>blue</i>	5070115
			1	5/0	75 cm	<i>blue</i>	5100115
			1.5	4/0	75 cm	<i>blue</i>	5150115
HR 17 mm			1	5/0	75 cm	<i>blue</i>	5100117
			2	3/0	75 cm	<i>blue</i>	5200117
HR 20 mm			1	5/0	75 cm	<i>blue</i>	5100120
			1.5	4/0	75 cm	<i>blue</i>	5150120
			2	3/0	75 cm	<i>blue</i>	5200120
HR 25 mm			3	2/0	75 cm	<i>blue</i>	5300120
			1	5/0	75 cm	<i>blue</i>	5100125
			1.5	4/0	75 cm	<i>blue</i>	5150125
			2	3/0	75 cm	<i>blue</i>	5200125
			3	2/0	75 cm	<i>blue</i>	5300125
HR 30 mm			3.5	0	75 cm	<i>blue</i>	5350125
			4	1	75 cm	<i>blue</i>	5400125
			2	3/0	75 cm	<i>blue</i>	5200130
			3	2/0	75 cm	<i>blue</i>	5300130
HR 35 mm			3.5	0	75 cm	<i>blue</i>	5350130
			4	1	75 cm	<i>blue</i>	5400130
			3	2/0	75 cm	<i>blue</i>	5300135
HR 40 mm			3.5	0	75 cm	<i>blue</i>	5350135
			4	1	75 cm	<i>blue</i>	5400135
			4	1	90 cm	<i>blue</i>	5400140
			5	2	75 cm	<i>blue</i>	5500140




Round bodied taper point


3/8 circle





			Metric	USP	Length		Order number
DR 16 mm			1.5	4/0	75 cm	<i>blue</i>	5150216





**Round bodied taper point
double armed**



		Metric	USP	Length		Order number
2 x HR 8 mm		0.4	8/0	60 cm	<i>blue</i>	5040108-2
		0.5	7/0	60 cm	<i>blue</i>	5050108-2
		0.7	6/0	60 cm	<i>blue</i>	5070108-2
2 x HR 9.3 mm		0.5	7/0	60 cm	<i>blue</i>	5050193-2
		0.7	6/0	60 cm	<i>blue</i>	5070193-2
2 x HR 10 mm		0.5	7/0	60 cm	<i>blue</i>	5050110-2
		0.7	6/0	60 cm	<i>blue</i>	5070110-2
		1	5/0	60 cm	<i>blue</i>	5100110-2
2 x HR 13 mm		0.5	7/0	60 cm	<i>blue</i>	5050113-2
		0.7	6/0	60 cm	<i>blue</i>	5070113-2
		1	5/0	75 cm	<i>blue</i>	5100113-2
2 x HR 16 mm		0.7	6/0	90 cm	<i>blue</i>	5070116-2
		1	5/0	90 cm	<i>blue</i>	5100116-2
		1.5	4/0	90 cm	<i>blue</i>	5150116-2
		2	3/0	90 cm	<i>blue</i>	5200116-2
2 x HR 20 mm		1	5/0	90 cm	<i>blue</i>	5100120-2
		1.5	4/0	90 cm	<i>blue</i>	5150120-2
		2	3/0	90 cm	<i>blue</i>	5200120-2
2 x HR 25 mm		1.5	4/0	90 cm	<i>blue</i>	5150125-2
		2	3/0	90 cm	<i>blue</i>	5200125-2
		3	2/0	90 cm	<i>blue</i>	5300125-2
2 x HR 30 mm		2	3/0	90 cm	<i>blue</i>	5200130-2
		3	2/0	90 cm	<i>blue</i>	5300130-2
2 x HR 35 mm		3	2/0	90 cm	<i>blue</i>	5300135-2
2 x HR 40 mm		3.5	0	90 cm	<i>blue</i>	5350140-2
		4	1	90 cm	<i>blue</i>	5400140-2

 Round bodied taper point double armed	Metric	USP	Length	Order number	
1/2 circle 	HR 15 mm	2	3/0	90 cm	<i>blue</i> 5200115-30
	HR 30 mm				

 Round bodied taper point double armed	Metric	USP	Length	Order number
2 x DR 6.5 mm	0.4	8/0	60 cm	<i>blue</i> 50402065-2-60
2 x DR 8 mm	0.4	8/0	45 cm	<i>blue</i> 5040208-2
	0.5	7/0	60 cm	<i>blue</i> 5050208-2
2 x DR 9.3 mm	0.5	7/0	60 cm	<i>blue</i> 5050293-2
	0.7	6/0	60 cm	<i>blue</i> 5070293-2
2 x DR 10 mm	0.7	6/0	60 cm	<i>blue</i> 5070210-2
	1	5/0	75 cm	<i>blue</i> 5100210-2
2 x DR 13 mm	0.5	7/0	60 cm	<i>blue</i> 5050213-2
	1	5/0	60 cm	<i>blue</i> 5100213-2
	0.7	6/0	60 cm	<i>blue</i> 5070213-2
3/8 circle 	1	5/0	90 cm	<i>blue</i> 5100218-2
	1.5	4/0	90 cm	<i>blue</i> 5150218-2
2 x DR 22 mm	1.5	4/0	90 cm	<i>blue</i> 5150222-2
2 x DR 26 mm	1.5	4/0	90 cm	<i>blue</i> 5150226-2
2 x DR 30 mm	3	2/0	90 cm	<i>blue</i> 5300230-2



Reverse cutting

Metric

USP



Length



Order number



3/8 circle













			Metric	USP	Length		Order number
DS 10 mm			0.7	6/0	45 cm	blue	5071510
DS 12 mm			0.5	7/0	45 cm	blue	5051512
			0.7	6/0	45 cm	blue	5071512
			1	5/0	45 cm	blue	5101512
			1.5	4/0	75 cm	blue	5151512
DS 12 mm Easy Pass			0.7	6/0	50 cm	blue	5071512P
			1	5/0	50 cm	blue	5101512P
DS 16 mm			0.7	6/0	45 cm	blue	5071516
			1	5/0	45 cm	blue	5101516
			1.5	4/0	45 cm	blue	5151516-45
			1.5	4/0	75 cm	blue	5151516
DS 16 mm Easy Pass			2	3/0	75 cm	blue	5201516
			0.7	6/0	50 cm	blue	5071516P
			1	5/0	50 cm	blue	5101516P
DS 16 mm black			1.5	4/0	50 cm	blue	5151516P
			1	5/0	45 cm	blue	5101516VB
			0.7	6/0	75 cm	blue	5071519
DS 19 mm			1	5/0	45 cm	blue	5101519-45
			1	5/0	75 cm	blue	5101519
			1.5	4/0	45 cm	blue	5151519-45
			1.5	4/0	75 cm	blue	5151519
			2	3/0	45 cm	blue	5201519-45
			2	3/0	75 cm	blue	5201519
			3	2/0	75 cm	blue	5301519
			3.5	0	75 cm	blue	5351519
DS 19 mm Easy Pass			1	5/0	50 cm	blue	5101519P
			1.5	4/0	50 cm	blue	5151519P
			2	3/0	50 cm	blue	5201519P
DS 19 mm black			1	5/0	75 cm	blue	5101519VB
			1.5	4/0	75 cm	blue	5151519VB
DS 24 mm			1.5	4/0	75 cm	blue	5151524
			2	3/0	45 cm	blue	5201524-45
			2	3/0	75 cm	blue	5201524
			3	2/0	45 cm	blue	5301524-45
			3	2/0	75 cm	blue	5301524
			3.5	0	75 cm	blue	5351524
	4	1	75 cm	blue	5401524		








 Reverse cutting		Metric	USP	Length	Order number
	DS 24 mm Easy Pass	1.5	4/0	50 cm	<i>blue</i> 5151524P
		2	3/0	50 cm	<i>blue</i> 5201524P
		3	2/0	50 cm	<i>blue</i> 5301524P
	DS 26 mm	2	3/0	45 cm	<i>blue</i> 5201526-45
		2	3/0	75 cm	<i>blue</i> 5201526
		3	2/0	45 cm	<i>blue</i> 5301526-45
	DS 30 mm	3	2/0	75 cm	<i>blue</i> 5301526
		2	3/0	75 cm	<i>blue</i> 5201530
		3	2/0	75 cm	<i>blue</i> 5301530
	DS 36 mm	3.5	0	75 cm	<i>blue</i> 5351530
		4	1	75 cm	<i>blue</i> 5401530
		2	3/0	75 cm	<i>blue</i> 5201536
3		2/0	75 cm	<i>blue</i> 5301536	
DS 40 mm	3.5	0	75 cm	<i>blue</i> 5351536	
	4	1	75 cm	<i>blue</i> 5401536	
	5	2	90 cm	<i>blue</i> 5501540	
DS 45 mm	3	2/0	75 cm	<i>blue</i> 5301540	
DS 45 mm	4	1	90 cm	<i>blue</i> 5401545	
DS 90 mm	4	1	90 cm	<i>blue</i> 5401590	

 Reverse cutting		Metric	USP	Length	Order number
	HS 22 mm	2	3/0	45 cm	<i>blue</i> 5201422
	HS 24 mm	3	2/0	75 cm	<i>blue</i> 5301424

 Cutting		Metric	USP	Length	Order number
	GS 60 mm	2	3/0	75 cm	<i>blue</i> 5202160
	GS 75 mm	3	2/0	75 cm	<i>blue</i> 5302160
	GS 75 mm	3	2/0	75 cm	<i>blue</i> 5302175

 Round bodied taper cutting point double armed		Metric	USP	Length	Order number
	2 x HRT 16 mm 	1	5/0	75 cm	<i>blue</i> 5100516-2
		1.5	4/0	90 cm	<i>blue</i> 5150516-2
	2 x HRT 22 mm 	1.5	4/0	90 cm	<i>blue</i> 5150522-2
		2 x HRT 25 mm 	1.5	4/0	90 cm
	2		3/0	90 cm	<i>blue</i> 5200525-2
	3		2/0	90 cm	<i>blue</i> 5300525-2
	2 x HRT 30 mm 	2	3/0	90 cm	<i>blue</i> 5200530-2
		3	2/0	90 cm	<i>blue</i> 5300530-2

 Round bodied taper cutting point		Metric	USP	Length	Order number
	HRT 10 mm 	0.7	6/0	75 cm	<i>blue</i> 5070510
	HRT 25 mm 	3	2/0	75 cm	<i>blue</i> 5300525

 Soft cutting point coronary double armed		Metric	USP	Length	Order number	
	2 x CC 6.4 mm 	0.4	8/0	60 cm	<i>blue</i> 5042606-2	
		0.4	8/0	60 cm	<i>blue</i> 5042608-2	
	2 x CC 8 mm 	0.5	7/0	75 cm	<i>blue</i> 5052608-2	
		0.7	6/0	75 cm	<i>blue</i> 5072608-2	
	2 x CC 9.3 mm 	0.5	7/0	60 cm	<i>blue</i> 5052693-2	
		0.7	6/0	60 cm	<i>blue</i> 5072693-2	
	2 x CC 10 m 	0.5	7/0	75 cm	<i>blue</i> 5052610-2	
		0.7	6/0	75 cm	<i>blue</i> 5072610-2	
	2 x CC 13 mm 	0.5	7/0	75 cm	<i>blue</i> 5052613-2	
		0.7	6/0	75 cm	<i>blue</i> 5072613-2	
			1	5/0	75 cm	<i>blue</i> 5102613-2



Soft cutting point coronary
black double armed

3/8 circle



		Metric	USP	Length		Order number	
2 x CC 8 mm black		0.4	8/0	60 cm	<i>blue</i>	5042608-2VB	
		0.5	7/0	75 cm	<i>blue</i>	5052608-2VB	
		0.7	6/0	75 cm	<i>blue</i>	5072608-2VB	
	2 x CC 9 mm black		0.4	8/0	60 cm	<i>blue</i>	5042609-2VB
			0.5	7/0	60 cm	<i>blue</i>	5052609-2VB
			0.7	6/0	75 cm	<i>blue</i>	5072613-2VB
2 x CC 13 mm black		1	5/0	75 cm	<i>blue</i>	5102613-2VB	